PORTUGUÊS 

Plenum Castable 2.0

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ENGLISH 

Plenum Castable 2.0

Instructions for Use

ESPAÑOL 

Plenum Castable 2.0

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Rev00

Plenum Castable 2.0

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Plenum Castable 2.0

Instructions for Use

Power, endurance,
and control.
No compromises.

Plenum Castable 2.0

Indications for use

Plenum Castable 2.0 is a light-curable resin designed for use in the production of castable components for casting and injection molding processes. This material is an alternative to casting wax.

Contraindications

This material is contraindicated for any purpose other than its intended use and when there is direct intraoral contact with resin that is not fully cured.

Device description

Plenum Castable 2.0 is a resin designed for printing castable parts, which are intended exclusively for professional dental use.

Manufacturing process

1. CAD File: File in STL format.

2. Prepare:

- 2.1** Make sure the printer you are using is compatible with Plenum Castable 2.0.
- 2.2** Set the recommended settings for the printer.
- 2.3** Sprintray Printers (Pro 95, Pro 55, and Pro 2): Use the Sprintray BR Model Beige print profile.

3. Shake:

- 3.1** Shake the bottle well before use (2-5 minutes).
- 3.2** Pour the resin into the tank until it reaches at least the minimum level mark.
- 3.3** Gently mix the resin in the tank before each print.

4. Printing:

- 4.1** Print the selected documents.

5. Remove:

- 5.1** Remove printed objects from the platform immediately after printing is complete.
- 5.2** Remove all supports using a cutter or cutting disc, if necessary.

6. Cleaning:

- 6.1** Isopropyl alcohol, 91% or higher.
- 6.2** Ultrasonic bath (Duration: 2 cycles, 3 minutes each).
- 6.3** ProWash S or ProWash/Dry: not recommended, as the wax content may cause the equipment to clog.

7. Drying:

- 7.1** Dry the printed devices with compressed air.
- 7.2** If compressed air is not available, let the devices dry naturally.

8. UV Post-Curing:

- 8.1** Full polymerization in the post-curing unit is essential to achieve the desired material properties.
- 8.2** 5 min in a 405 nm UV post-cure unit.
- 8.3** Power of 90W or higher.
- 8.4** ProCure 1: Castable 2.
- 8.5** ProCure 2: Castable 2.
- 8.6** NanoCure: Castable 2.

9. Burning:

- 9.1** Insert the lining into the ring.
- 9.2** Preheat the oven to 800°C.
- 9.3** After the lining has set (approx. 10 to 20 minutes), insert the ring at a slight angle with the opening facing down; if necessary, use a small piece of lining to serve as a support.
- 9.4** Leave the ring at 800°C for 20 minutes.
- 9.5** Set the temperature to 930°C and let it sit for another 50 minutes, starting the timer from the moment you set the new temperature.
- 9.6** Proceed to the injection process, paying close attention to the temperatures recommended by the alloy or disilicate manufacturer for this step.

Warnings and precautions

Plenum Castable 2.0 is not a biocompatible material in its processed and cured form. In its uncured form, it is classified as a sensitizing material. When washing with a solvent, do so in a well-ventilated area while wearing appropriate protective equipment.

- **Skin contact:** May cause skin irritation. If unprocessed resin comes into contact with the skin, wash thoroughly with soap and water. May cause an allergic skin reaction. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical attention.
- **Inhalation:** High concentrations of vapor may cause headaches, eye irritation, and/or respiratory irritation. If exposed to high concentrations of vapor or mist, move to fresh air. Use oxygen or artificial respiration as needed.
- **Eye contact:** Thoroughly wash the affected area with soap and water.
- **Ingestion:** Contact a doctor immediately.

Storage

- **Reusing Material:** The resin remaining in the resin tank can be reused. You can use a filter to ensure that the resin is free of cured particles to prevent print failures. The remaining material in the tank can be poured back into the resin bottle after filtering. This process can be repeated until the material in the bottle is completely used up. Please note that when reusing the resin, it must be filtered and poured back into the same bottle.
- Store Plenum Castable 2.0 at temperatures between 2°C and 28°C (35.6°F and 82.4°F) and keep it out of direct sunlight.
- Keep the bottle closed and/or the tank cap securely fastened when not in use.

- Allow the product to cure completely before disposal.
- Do not use Plenum Castable 2.0 after the expiration date printed on the bottle.

Do not use expired resin; print stability may be compromised if expired photoinitiators are not activated properly.

Contact Information

Manufacturer: ASTRO SCIENCE DO BRASIL S.A.

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Fabricante: ASTRO SCIENCE DO BRASIL S.A.
Ventas y Distribución: M3 Health Industria e Comercio de Produtos Médicos, Odontológicos e Correlatos S.A