PORTUGUÊS 

Plenum Crown

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ENGLISH 

Plenum Crown

Instructions for Use

ESPAÑOL 

Plenum Crown

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Rev00

Plenum Crown

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Plenum Crown

Instructions for Use

Power, endurance,
and control.
No compromises.

Plenum Crown

Indications for use

Plenum Crown is a light-curing resin designed for use in the fabrication of single crowns, veneers, bridges, inlays, and onlays.

Contraindications

This material is contraindicated for any purpose other than its intended use and when there is direct intraoral contact with resin that is not fully cured.

Device description

Plenum Crown is a Class II resin intended for printing dental restorations with a long service life, which are to be used exclusively for professional dental work.

Manufacturing process

1. CAD File: A file in STL format with the following characteristics:

- 1.1** Crown walls: ≥ 1.0 mm
- 1.2** Crown margins: ≥ 1.0 mm
- 1.3** Lens thickness: ≥ 0.5 mm

2. Prepare:

- 2.1** Make sure the printer you are using is compatible with Plenum Crown.
- 2.2** Set the recommended settings for the printer.
- 2.3** Sprinray Printers (Pro 95, Pro 55, and Pro 2): Use the Sprinray Crown (Bego) print profile.

3. Shake:

- 3.1** Shake the bottle well before use (10–15 minutes). If necessary, place the bottle in a container of warm water to facilitate complete homogenization of the ceramic filler.
- 3.2** Pour the resin into the tank until it reaches at least the minimum level mark.
- 3.3** Gently mix the resin in the tank before each print.

4. Printing:

- 4.1** Print the selected documents.

5. Remove:

- 5.1** Remove printed objects from the platform immediately after printing is complete.
- 5.2** Remove all supports using a cutter or cutting disc, if necessary.

6. Cleaning:

- 6.1** Isopropyl alcohol, 91% or higher.
- 6.2** Wash the device: submerge it in a small container and scrub it with a brush.
- 6.3** Pro Wash S or Pro Wash/Dry: We recommend using cycles of 3 to 4 minutes at a time.

7. Drying:

- 7.1** Dry the printed devices with compressed air.
- 7.2** If compressed air is not available, let the devices dry naturally.

8. UV Post-Curing:

- 8.1** Full polymerization in the post-curing unit is essential to achieve the desired material properties.
- 8.2** 60 minutes in a 405-nm UV post-cure unit.
- 8.3** Power of 90W or higher.
- 8.4** ProCure 1: Bego (Crown).
- 8.5** ProCure 2: Sprintray Crown (Bego).
- 8.6** NanoCure: Sprintary Crown (Bego).

Warnings and precautions

Plenum Crown is non-toxic in its processed and cured form and is classified as a biocompatible material. In its uncured form, it is classified as a sensitizing material.

When washing with a solvent, do so in a well-ventilated area while wearing appropriate protective equipment.

- **Skin contact:** May cause skin irritation. If unprocessed resin comes into contact with the skin, wash thoroughly with soap and water. May cause an allergic skin reaction. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical attention.
- **Inhalation:** High concentrations of vapor may cause headaches, eye irritation, and/or respiratory irritation. If exposed to high concentrations of vapor or mist, move to fresh air. Use oxygen or artificial respiration as needed.
- **Eye contact:** Thoroughly wash the affected area with soap and water.
- **Ingestion:** Contact a doctor immediately.

Storage

- **Reusing Material:** The resin remaining in the resin tank can be reused. You can use a filter to ensure that the resin is free of cured particles to prevent print failures. The remaining material in the tank can be poured back into the resin bottle after filtering. This process can be repeated until the material in the bottle is completely used up. Please note that when reusing the resin, it must be filtered and poured back into the same bottle.
- Store Plenum Crown between 2°C and 28°C (35.6°F and 82.4°F) and keep it out of direct sunlight.
- Keep the bottle closed and/or the tank cap securely fastened when not in use.
- Allow to cure completely before disposal.
- Do not use Plenum Crown after the expiration date printed on the bottle.

Do not use expired resin; print stability may be compromised if expired photoinitiators are not activated properly.

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Fabricante / Titular del registro: ASTRO SCIENCE DO BRASIL S.A.
Ventas y Distribución: M3 Health Industria e Comercio de Produtos Médicos, Odontológicos e Correlatos S.A