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HEALING
ABUTMENT Ti
English

CICATRIZADOR Ti
Español

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HEALING ABUTMENT Ti

INSTRUCTIONS FOR USE

TECHNICAL NAME: Dental Implant Components

TRADE NAME: Healing abutment Ti

MODEL: Detailed product description see table 1

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E-mail: relacionamento@plenum.bio

Website: www.plenum.bio

ATTENTION: Read all Instructions for Use carefully before use.

Comply with all the indications, warnings and precautions mentioned in these Instructions for Use. Failure to observe these points could lead to complications. Product for Medical Use. Sterile product. Sterilized by Gamma Irradiation. Shelf life: See product labeling. Manufacturer recommends Single Use. Do not resterilize and do not reuse.

ANVISA REGISTRATION No.: 81684340022

Code, Batch Number, Date of Manufacture and Sterilization: See product labeling.

Rev.01

SUMMARY

1. DETAILED DESCRIPTION OF THE MEDICAL PRODUCT	3
2. PRODUCT COMPOSITION	4
3. FORM OF PRESENTATION OF THE MEDICAL PRODUCT	4
4. INDICATION, PURPOSE OR USE FOR WHICH THE MEDICAL PRODUCT IS INTENDED	5
5. FUNDAMENTALS OF THE PRODUCT’S OPERATION AND INTENDED PERFORMANCE	5
6. INSTRUCTIONS FOR USING THE PRODUCT	5
7. PRECAUTIONS	5
8. WARNINGS	6
9. CONTRAINDICATIONS	6
10. MAGNETIC RESONANCE IMAGING (MRI) - INFORMATION	7
11. POSSIBLE ADVERSE EFFECTS	7
12. SYMBOLS	7
13. AUXILIARY INSTRUMENTS	8
14. STERILIZATION	8
15. SHELF LIFE OF THE PRODUCT	8
16. IDENTIFICATION AND TRACEABILITY	8
17. CONDITIONS STORAGE STORAGE CONDITIONS, TRANSPORTATION, CONSERVATION AND HANDLING	8
18. DISPOSAL OF THE PRODUCT	8
19. CUSTOMER SERVICE	9
20. LIMITED WARRANTY	9

1. DETAILED DESCRIPTION OF THE MEDICAL PRODUCT

The Healing abutment Ti is a temporary medical device to be installed exclusively on Plenum Implants in order to guide the proper healing or remodeling of the gingival tissue for the placement of a dental prosthesis, after implant placement or during reopening surgery (second surgical phase), shaping the gingival tissue in order to maintain a space for late dental prosthesis placement procedures. The component height is indicated according to the patient’s interocclusal space.

The fitting and fixation between the Healing Cap and the implant is done through a thread and Morse taper type connection. The upper part of the Healing Cap has a fitting for the 1.2mm hexagonal key, a connection for installing the product on the dental implant.

The products that are part of the Healing abutment Ti family are available in 5 (five) prosthetic adaptations: Slim (SL) represented by the color pink, Regular (RE) represented by the color yellow, Unifit (UF) represented by the color gray,

Morsefit (MF) represented by the color blue and Short (SH) represented by the color green. They must be used in conjunction with the Plenum Dental Implant and its respective prosthetic adapter.

Each commercial code of the products that make up the Healing abutment Ti family is composed of the following parts represented in Figure 1 and as per Table 1.

The Healing abutment Ti in the models described in Table 1 should only be used with Plenum brand Dental Implants and Instruments (these Implants and Instruments are not part of these instructions for use and must be purchased separately).


FIGURE 1: Schematic model of Ti scarers.

TABLE 1: Illustrative images, codes, models, prosthetic platforms and dimensions of Healing abutment Tis

Imagem ilustrativa	Código	Nome / Modelo comercial	Composição do Modelo Comercial	Diâmetro (mm)	Comprimento (mm)	Interface protética
	018-21-C00-35-02	HEALING ABUTMENT TI SL Ø3.5 X 2MM	HEALING ABUTMENT TI SL Ø3.5 X 2MM	3.5	2.00	SL
	018-21-C00-35-03	HEALING ABUTMENT TI SL Ø3.5 X 3MM	HEALING ABUTMENT TI SL Ø3.5 X 3MM	3.5	3.00	SL
	018-21-C00-35-04	HEALING ABUTMENT TI SL Ø3.5 X 4MM	HEALING ABUTMENT TI SL Ø3.5 X 4MM	3.5	4.00	SL
	018-21-C00-35-05	HEALING ABUTMENT TI SL Ø3.5 X 5MM	HEALING ABUTMENT TI SL Ø3.5 X 5MM	3.5	5.00	SL
	018-21-C00-35-06	HEALING ABUTMENT TI SL Ø3.5 X 6MM	HEALING ABUTMENT TI SL Ø3.5 X 6MM	3.5	6.00	SL
	018-22-C00-35-02	HEALING ABUTMENT TI RE Ø3.5 X 2MM	HEALING ABUTMENT TI RE Ø3.5 X 2MM	3.5	2.05	RE
	018-22-C00-35-03	HEALING ABUTMENT TI RE Ø3.5 X 3MM	HEALING ABUTMENT TI RE Ø3.5 X 3MM	3.5	3.05	RE
	018-22-C00-35-04	HEALING ABUTMENT TI RE Ø3.5 X 4MM	HEALING ABUTMENT TI RE Ø3.5 X 4MM	3.5	4.05	RE
	018-22-C00-35-05	HEALING ABUTMENT TI RE Ø3.5 X 5MM	HEALING ABUTMENT TI RE Ø3.5 X 5MM	3.5	5.05	RE
	018-22-C00-35-06	HEALING ABUTMENT TI RE Ø3.5 X 6MM	HEALING ABUTMENT TI RE Ø3.5 X 6MM	3.5	6.05	RE
	018-22-C00-45-02	HEALING ABUTMENT TI RE Ø4.5 X2MM	HEALING ABUTMENT TI RE Ø4.5 X 2MM	4.5	2.05	RE
	018-22-C00-45-03	HEALING ABUTMENT TI RE Ø4.5 X3MM	HEALING ABUTMENT TI RE Ø4.5 X 3MM	4.5	3.05	RE
	018-22-C00-45-04	HEALING ABUTMENT TI RE Ø4.5 X4MM	HEALING ABUTMENT TI RE Ø4.5 X 4MM	4.5	4.05	RE
	018-22-C00-45-05	HEALING ABUTMENT TI RE Ø4.5 X5MM	HEALING ABUTMENT TI RE Ø4.5 X 5MM	4.5	5.05	RE
	018-22-C00-45-06	HEALING ABUTMENT TI RE Ø4.5 X 6MM	HEALING ABUTMENT TI RE Ø4.5 X 6MM	4.5	6.05	RE
	018-22-C00-55-02	HEALING ABUTMENT TI RE Ø5.5 X 2MM	HEALING ABUTMENT TI RE Ø5.5 X 2MM	5.5	2.05	RE
	018-22-C00-55-03	HEALING ABUTMENT TI RE Ø5.5 X 3MM	HEALING ABUTMENT TI RE Ø5.5 X 3MM	5.5	3.05	RE
	018-22-C00-55-04	HEALING ABUTMENT TI RE Ø5.5 X 4MM	HEALING ABUTMENT TI RE Ø5.5 X 4MM	5.5	4.05	RE
	018-22-C00-55-05	HEALING ABUTMENT TI RE Ø5.5 X 5MM	HEALING ABUTMENT TI RE Ø5.5 X 5MM	5.5	5.05	RE
	018-22-C00-55-06	HEALING ABUTMENT TI RE Ø5.5 X 6MM	HEALING ABUTMENT TI RE Ø5.5 X 6MM	5.5	6.05	RE
	018-24-C00-35-02	HEALING ABUTMENT TI UF Ø3.5 X 2MM	HEALING ABUTMENT TI UF Ø3.5 X 2MM	3.5	2.00	UF
	018-24-C00-35-03	HEALING ABUTMENT TI UF Ø3.5 X 3MM	HEALING ABUTMENT TI UF Ø3.5 X 3MM	3.5	3.00	UF
	018-24-C00-35-04	HEALING ABUTMENT TI UF Ø3.5 X 4MM	HEALING ABUTMENT TI UF Ø3.5 X 4MM	3.5	4.00	UF
	018-24-C00-35-05	HEALING ABUTMENT TI UF Ø3.5 X 5MM	HEALING ABUTMENT TI UF Ø3.5 X 5MM	3.5	5.00	UF
	018-24-C00-35-06	HEALING ABUTMENT TI UF Ø3.5 X 6MM	HEALING ABUTMENT TI UF Ø3.5 X 6MM	3.5	6.00	UF
	018-24-C00-45-02	HEALING ABUTMENT TI UF Ø4.5 X 2MM	HEALING ABUTMENT TI UF Ø4.5 X 2MM	4.5	2.00	UF
	018-24-C00-45-03	HEALING ABUTMENT TI UF Ø4.5 X 3MM	HEALING ABUTMENT TI UF Ø4.5 X 3MM	4.5	3.00	UF
	018-24-C00-45-04	HEALING ABUTMENT TI UF Ø4.5 X 4MM	HEALING ABUTMENT TI UF Ø4.5 X 4MM	4.5	4.00	UF
	018-24-C00-45-05	HEALING ABUTMENT TI UF Ø4.5 X 5MM	HEALING ABUTMENT TI UF Ø4.5 X 5MM	4.5	5.00	UF
	018-24-C00-45-06	HEALING ABUTMENT TI UF Ø4.5 X 6MM	HEALING ABUTMENT TI UF Ø4.5 X 6MM	4.5	6.00	UF
	018-24-C00-55-02	HEALING ABUTMENT TI UF Ø5.5 X 2MM	HEALING ABUTMENT TI UF Ø5.5 X 2MM	5.5	2.00	UF
	018-24-C00-55-03	HEALING ABUTMENT TI UF Ø5.5 X 3MM	HEALING ABUTMENT TI UF Ø5.5 X 3MM	5.5	3.00	UF

	018-24-C00-55-04	HEALING ABUTMENT TI UF Ø5.5 X 4MM	HEALING ABUTMENT TI UF Ø5.5 X 4MM	5.5	4.00	UF
	018-24-C00-55-05	HEALING ABUTMENT TI UF Ø5.5 X 5MM	HEALING ABUTMENT TI UF Ø5.5 X 5MM	5.5	5.00	UF
	018-24-C00-55-06	HEALING ABUTMENT TI UF Ø5.5 X 6MM	HEALING ABUTMENT TI UF Ø5.5 X 6MM	5.5	6.00	UF
	018-25-C00-35-02	HEALING ABUTMENT TI MF Ø3.5 X 2MM	HEALING ABUTMENT TI MF Ø3.5 X 2MM	3.5	2.00	MF
	018-25-C00-35-03	HEALING ABUTMENT TI MF Ø3.5 X 3MM	HEALING ABUTMENT TI MF Ø3.5 X 3MM	3.5	3.00	MF
	018-25-C00-35-04	HEALING ABUTMENT TI MF Ø3.5 X 4MM	HEALING ABUTMENT TI MF Ø3.5 X 4MM	3.5	4.00	MF
	018-25-C00-35-05	HEALING ABUTMENT TI MF Ø3.5 X 5MM	HEALING ABUTMENT TI MF Ø3.5 X 5MM	3.5	5.00	MF
	018-25-C00-35-06	HEALING ABUTMENT TI MF Ø3.5 X 6MM	HEALING ABUTMENT TI MF Ø3.5 X 6MM	3.5	6.00	MF
	018-25-C00-45-02	HEALING ABUTMENT TI MF Ø4.5 X 2MM	HEALING ABUTMENT TI MF Ø4.5 X 2MM	4.5	2.00	MF
	018-25-C00-45-03	HEALING ABUTMENT TI MF Ø4.5 X 3MM	HEALING ABUTMENT TI MF Ø4.5 X 3MM	4.5	3.00	MF
	018-25-C00-45-04	HEALING ABUTMENT TI MF Ø4.5 X 4MM	HEALING ABUTMENT TI MF Ø4.5 X 4MM	4.5	4.00	MF
	018-25-C00-45-05	HEALING ABUTMENT TI MF Ø4.5 X 5MM	HEALING ABUTMENT TI MF Ø4.5 X 5MM	4.5	5.00	MF
	018-25-C00-45-06	HEALING ABUTMENT TI MF Ø4.5 X 6MM	HEALING ABUTMENT TI MF Ø4.5 X 6MM	4.5	6.00	MF
	018-25-C00-55-02	HEALING ABUTMENT TI MF Ø5.5 X 2MM	HEALING ABUTMENT TI MF Ø5.5 X 2MM	5.5	2.00	MF
	018-25-C00-55-03	HEALING ABUTMENT TI MF Ø5.5 X 3MM	HEALING ABUTMENT TI MF Ø5.5 X 3MM	5.5	3.00	MF
	018-25-C00-55-04	HEALING ABUTMENT TI MF Ø5.5 X 4MM	HEALING ABUTMENT TI MF Ø5.5 X 4MM	5.5	4.00	MF
	018-25-C00-55-05	HEALING ABUTMENT TI MF Ø5.5 X 5MM	HEALING ABUTMENT TI MF Ø5.5 X 5MM	5.5	5.00	MF
	018-25-C00-55-06	HEALING ABUTMENT TI MF Ø5.5 X 6MM	HEALING ABUTMENT TI MF Ø5.5 X 6MM	5.5	6.00	MF
	018-23-C00-45-02	HEALING ABUTMENT TI SH Ø4.5 X 2MM	HEALING ABUTMENT TI SH Ø4.5 X 2MM	4.5	2.05	SH
	018-23-C00-45-03	HEALING ABUTMENT TI SH Ø4.5 X 3MM	HEALING ABUTMENT TI SH Ø4.5 X 3MM	4.5	3.05	SH
	018-23-C00-45-04	HEALING ABUTMENT TI SH Ø4.5 X 4MM	HEALING ABUTMENT TI SH Ø4.5 X 4MM	4.5	4.05	SH
	018-23-C00-45-05	HEALING ABUTMENT TI SH Ø4.5 X 5MM	HEALING ABUTMENT TI SH Ø4.5 X 5MM	4.5	5.05	SH
	018-23-C00-45-06	HEALING ABUTMENT TI SH Ø4.5 X 6MM	HEALING ABUTMENT TI SH Ø4.5 X 6MM	4.5	6.05	SH
	018-23-C00-55-02	HEALING ABUTMENT TI SH Ø5.5 X 2MM	HEALING ABUTMENT TI SH Ø5.5 X 2MM	5.5	2.05	SH
	018-23-C00-55-03	HEALING ABUTMENT TI SH Ø5.5 X 3MM	HEALING ABUTMENT TI SH Ø5.5 X 3MM	5.5	3.05	SH
	018-23-C00-55-04	HEALING ABUTMENT TI SH Ø5.5 X 4MM	HEALING ABUTMENT TI SH Ø5.5 X 4MM	5.5	4.05	SH
	018-23-C00-55-05	HEALING ABUTMENT TI SH Ø5.5 X 5MM	HEALING ABUTMENT TI SH Ø5.5 X 5MM	5.5	5.05	SH
	018-23-C00-55-06	HEALING ABUTMENT TI SH Ø5.5 X 6MM	HEALING ABUTMENT TI SH Ø5.5 X 6MM	5.5	6.05	SH

NOTE: The component colors represent the prosthetic interface: PINK for the SLIM prosthetic connection platform interface; YELLOW for the REGULAR prosthetic connection platform interface; GRAY for the UNIFIT prosthetic connection platform interface; BLUE for the Morsefit prosthetic connection platform interface; and GREEN for the SHORT prosthetic connection platform interface.

2. PRODUCT COMPOSITION

The models in the Healing abutment Ti family are made from Pure Grade 4 Titanium, in accordance with ASTM F67 and ISO 5832-2 specifications.


3. FORM OF PRESENTATION OF THE MEDICAL PRODUCT

The Healing abutment Ti in the models described in Tables 1 are supplied unit-packed, in sterile condition, and presented as follows:

- Primary packaging: polyethylene terephthalate (PET) blister and heat-sealed with surgical grade paper; containing 01 unit of the product, sterilized using Gamma Irradiation;
- Secondary packaging: rigid carton containing the primary packaging and 03 identification and traceability labels;

Below are illustrations of the primary and secondary packaging of a model of Healing abutment Ti, as it is sold.

IMAGE ILLUSTRATION OF PACKAGING PRIMARY

ILLUSTRATIVE IMAGE OF SECONDARY PACKAGING

ACCESSORIES:

The Healing abutment Ti has no accessories.

4. INDICATION, PURPOSE OR USE FOR WHICH THE MEDICAL PRODUCT IS INTENDED

The Healing abutment Ti is indicated and used during the healing of the gum tissue after the installation of a Plenum Dental Implant or at the time of reopening to remove the implant cover, in the preparation and installation of a dental prosthesis. These products in their various models must be installed in an exclusive connection with Plenum Dental Implants, with a view to the proper placement of a dental prosthesis

Note: The Healing abutment Ti should only be used with Plenum Dental Implants and Instruments (products that are not part of these instructions for use and must be purchased separately and are subject to specific ANVISA Registration). Implants and instruments from other brands may cause dimensional incompatibility, improper fitting, micro-movements and increase the risk of corrosion (electrolytic or galvanic) and fracture of these products, compromising the outcome of the surgical and prosthetic procedure.

5. FUNDAMENTALS OF THE PRODUCT’S OPERATION AND INTENDED PERFORMANCE

The Healing abutment Ti in its various models has the function of guiding, maintaining and conditioning the proper healing of the peri-implant gingival tissue after implant installation, in the reopening surgery (second surgical phase), shaping the gingival tissue in order to maintain a space for procedures to place the late dental prosthesis. The height of the component body is indicated according to the patient’s interocclusal space.

6. INSTRUCTIONS FOR USING THE PRODUCT

- Read all instructions for use, recommendations and guidelines carefully before using Healing abutment Ti;
- Before rehabilitating the patient, the dental surgeon must carry out rigorous planning using clinical and imaging tests (radiographic and tomographic) to determine whether the implants installed are in a position to receive prosthetic components;
- Before starting clinical and/or surgical procedures, the dental surgeon must subject the patient to rigorous oral antisepsis to prevent the component from coming into contact with contaminating substances;
- Before using Healing abutment Ti, the dental surgeon must rigorously check that the product has the same prosthetic interface as the implant/prosthetic component being installed, according to the table;
- Before starting the clinical procedure, the dental surgeon must ensure that the prosthetic components and auxiliary instruments are intact and complete;
- The clinical/surgical procedure, the model and size of the component, the technique and sequence of installation of the Healing abutment Ti must comply

with established implant dentistry techniques;

- In order to use Healing abutment Ti in a one-stage procedure (esthetics and/or immediate loading), the dental surgeon must respect the principles of occlusion established in implant dentistry and use compatible restorative materials so that the procedure has satisfactory results;
- For the use of Healing abutment Ti in two-stage procedures, prior conditioning of the peri-implant tissues with this component is recommended;
- Only use instruments manufactured by Plenum for the use of Healing abutment Ti, which are considered auxiliary instruments, are not part of these instructions for use and must be purchased separately (see Auxiliary Instruments);
- The product must be handled with care, especially during the process of attaching the component to the dental implant. The component should rest passively on the implant platform, avoiding a forced fit, the use of excessive force or oblique insertion;
- The use of imaging scans is highly recommended to certify the correct seating of the prosthetic component on the implant platform;
- Special care must be taken with the hexagonal connection on some components and the key used to insert them, as the two are not properly fitted elements can cause damage to this connection, compromising the integrity of the product and proper fixation to the implant;
- Care must be taken with the thread surface and the Morse Cone connection on the prosthetic component. Both surfaces must be free of particles and solid materials that could prevent a perfect fit and fixation with the implant and dental prosthesis, thus preventing proper fixation of these elements;
- The dental surgeon must carefully evaluate the insertion and removal of a prosthetic component more than once in the same patient, as this can cause deterioration of the threads, the Morse Cone and the component connection;
- Check passivity and carry out occlusal and interproximal adjustment after installing the prosthesis, avoiding compromising the implant/prosthesis set.

7. PRECAUTIONS

The dental surgeon must not use this product before carefully reading all the instructions for use and the information on the packaging. Failure to observe these points could lead to complications for the dental surgeon and especially for the patient;

- Cicatrizador Ti should be discarded after the gum tissue has healed.
- The main function of the Ti Cicatrizadores family is to be used in tissue preparation and the installation of a dental prosthesis. The dental surgeon is responsible for choosing the product models with the most appropriate anatomical and dimensional characteristics, after individual clinical analysis of each case;
- Immediate and late post-operative care and the correct oral hygiene technique are extremely important. The patient must be instructed on the risks that exist for the healthy maintenance of the implants and components; the patient’s ability and willingness to follow the dental surgeon’s instructions is one of the most important aspects in a dental implant and component placement procedure;
- The patient should be warned that failure to follow the post-operative instructions could lead to failure and early loss of the component and implant, requiring further surgery to remove the component or component/implant assembly;
- For the safety and effectiveness of the treatment, only prosthetic components and implants manufactured by Plenum should be used. Components and/

or implants manufactured by other companies may cause dimensional incompatibility, inadequate fitting, micro-movements, and increase the risk of fracture of these products, compromising the outcome of the surgical and prosthetic procedure;

- Special attention should be paid to the insertion keys before inserting the Healing abutment Ti. Damage to the tip of the keys can prevent the correct installation of the component, compromising the good prosthetic result. Therefore, do not use this product if any specific instruments have been damaged;

- The wrong choice of component to be used, as well as errors in indication, handling and installation, can cause damage to the product's structure, which can contribute to failure due to fatigue, fracture or loosening;

- The use of excessive torque during the insertion of a prosthetic component can lead to product failure. Do not exceed 10 N.cm of insertion torque for the fixing screw;

- If a prosthetic component and/or implant is overloaded beyond its functional capacity, excessive bone loss, implant failure, breakage of products and other elements of the prosthetic restoration can occur;

- Do not force the insertion of a prosthetic component beyond the depth established by the thread limit, as this may cause deformation of the hexagonal connection of the component and/or the insertion key and fracture of the product;

- Incorrect handling of small components inside the patient's mouth may result in the risk of aspiration and/or swallowing;

- The dental surgeon must rigorously assess the initial mechanical stability of the prosthetic component. If optimal and safe mechanical stability has not been achieved, the dental surgeon should evaluate and consider its possible removal and/or replacement;

- The patient must be informed of the risks involved in applying excessive mechanical load or trauma affecting the operated region, which could result in failure, loss of mechanical stability, loosening, breakage, fracture or loss of the dental implant and its components;

- Elderly patients, those with mental health problems, drug addicts, or those with any other difficulty in following post-operative recommendations may represent a greater risk of failure of this product, as these patients may ignore the instructions and restrictions that come with this procedure;

- Smokers, tobacco users, drug users, alcoholics, diabetics or malnourished patients should be warned about the increased incidence of lack of osseointegration, with consequent failure of implants and components; in addition, for this group of patients the rate of failure of the healing of the gingival tissue around the component body is considerable and may also be inefficient, compromising the installation of the dental prosthesis;

- If there are any visible deformations, scratches, cracks or dents on the component, it must be discarded;

- The patient must be informed of all the limitations, risks, complications, warnings and possible adverse effects of the dental implant and prosthetic component installation procedure;

- The dental surgeon must clearly guide the patient, who must understand the importance of post-operative follow-up, in which the stability, integrity of the gingival tissues and positioning of the prosthetic component and dental implant used must be checked, as well as the osseointegration and biological fixation of the implant at the appropriate time;

- In order to protect the patient and their dental surgeon, it is important to ensure the correct and complete traceability of the component used, by writing down the product code and batch in the patient's protocol, or by attaching the identification and traceability labels that come with the product to the patient's protocol;

- The success of a component placement procedure is related to the correct indication, selection, strength of the implant and component, correct surgical technique, as well as the quality of osseointegration, which can be affected by numerous biological, clinical and biomechanical factors, which are extrinsic to the quality of the implants and components. These factors can limit the useful life of the products and consequently cause them to fail and be lost prematurely. Therefore, faithful compliance with all the indications for use, contraindications, precautions, warnings and other observations contained in these instructions for use are essential for the success of the procedure.

8. WARNINGS

- Single-use product. Do not resterilize or reuse.

- Do not use if the packaging is open or damaged or if the expiration date has passed.

9. CONTRAINDICATIONS

- Healing abutment Ti are always used in conjunction with Plenum Dental Implants (supplied separately). Therefore, the dental surgeon should prudently assess the patient's condition prior to surgery and use of the components, and check for any contraindications listed below. Contraindications may include, but are not limited to:

- An implant that is poorly positioned and/or has insufficient osseointegration to properly support or anchor a dental prosthesis and its components is contraindicated for use with this product;

- Patients with poor bone quality due to osteoporosis or osteopenia, which can lead to early failure of the dental implant and component due to lack of osseointegration;

- Patients showing signs of allergy or hypersensitivity to Titanium Grade 4.

- Insufficient interocclusal or mesio-distal space for the insertion of any model of Healing abutment Ti;

- Periodontal inflammation and/or infection;

- Recent history of systemic or localized infection;

- Parafunctional habits;

- Patient with insufficient intermandibular space;

- Untreatable occlusion and joint problems;

- Pregnancy and lactation;

- Use in patients with immunological alterations and systemic diseases that result in poor healing or lack of integrity of the existing tissue at the implant site and gingival tissue;

- An alcoholic or drug-abusing patient who cannot follow the post-operative instructions correctly;

- Patient with mental illness or unable to follow the recommendations and care of the post-operative phase;

- Patient unable to maintain good oral hygiene.

10. MAGNETIC RESONANCE IMAGING (MRI) - INFORMATION

The product contains metal parts that may be affected in a magnetic resonance imaging (MRI) environment.

Non-clinical tests and MRI simulations have demonstrated that Plenum products made of metallic material can be used conditionally in MRI (MR Conditional).

A patient with this device can be safely scanned in a magnetic resonance imaging (MRI) system under the following conditions:

Static Magnetic Field: ≤ 3.0 T

Maximum Spatial Gradient Magnetic Field: 4,000 gauss/cm (40 T/m)

Maximum Indication for MRI System - Whole-body human exposure Specific Absorption Rate (SAR): 2.3 W/kg (15 minutes of scanning, normal operating mode);

Under the scanning conditions defined above, the medical device is expected to produce a maximum temperature increase of 0.8 °C after 15 minutes of continuous scanning.

In non-clinical tests, an image artifact originating from the medical device extends approximately 11.3 mm from the device when displayed with a gradient echo pulse sequence in a 3T MRI system.

11. POSSIBLE ADVERSE EFFECTS

No adverse effects are expected as long as the product is used according to its instructions. However, positive results may not be achieved in all surgical and prosthetic procedures. Every procedure presents risks and the possibility of complications, such as infections, allergic drug reactions, among others. Possible adverse effects associated with the use of prosthetic components may include, but are not limited to:

- Pain or discomfort due to the presence of the prosthetic component and/or dental implant, or as a result of the surgical procedure;
- Lack of osseointegration, resorption of bone tissue, or delayed bone healing, which can result in failure and breakage of the prosthetic component and/or dental implant;
- Fracture, displacement or early or late failure of the component and/or implant due to inadequate installation technique;
- Failure, displacement or fracture of the component or implant due to excessive loads applied to the products and/or inadequate bone fixation;
- Loss of the component or implant due to parafunctional habits or poor oral hygiene;
- Allergic reaction to the raw material of the Ti, Titanium Healing Cap Grade 4;
- Inefficient healing process of the gingival tissue due to an imbalance in the normal forces/tensions on the implant/component interface or the appearance of systemic disease(s) after the installation of the component, which the patient did not report to the dental surgeon;

- Bone necrosis, insufficient vascularization, bone resorption and insufficient bone formation, which can cause loss of fixation of the dental implant and consequently of the component; these effects are correlated to the patient’s systemic factors or inadequate surgical procedure without compliance with the instructions for use described for this family of components;

- The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient must be informed by the dental surgeon of all these complications, and of the risks arising from implant surgery.

12. SYMBOLS

The table describes the symbols that may be printed on the product labeling (primary and secondary packaging). Please check the physical labeling to see which symbols apply to the product.

	Manufacturer
	Date and country of manufacture
	Validity
	Lot
	Catalog number
	Sterilized with irradiation
	Do not resterilize
	No use if a packaging is damaged
	Unique sterile barrier system
	Protect from sunlight
	Keep dry
	Do not reuse
	Consult the instructions for use
	Unique device identifier
	Raw materials
	CAUTION: Federal law restricts this device to sale by or on the order of a doctor or dentist. (Notification required by the FDA for marketing in the United States of America)
	No reprocessing (Brazilian requirement)
	Conditional use in Magnetic Resonance Imaging

NOTA: Os símbolos gráficos relacionados à rotulagem atendem aos requisitos estabelecidos na norma ISO 15223-1.

13. AUXILIARY INSTRUMENTS

- Auxiliary instruments are all surgical and prosthetic instruments developed exclusively for the installation of the Healing abutment Ti;
- Plenum instruments are made from stainless steel (AISI 420-10), whose raw material provides the products with high resistance and durability, in accordance with the requirements specified by standards NBR 13911; ASTM F899 and ISO 5832-1;
- It is essential to use instruments manufactured by Plenum to guarantee the safe and effective installation of implants in accordance with the surgical protocol recommended by the manufacturer;
- It is essential that the dental surgeon uses only Plenum-made insertion driver to avoid damaging the component connection and the implant driver;
- For the proper installation of implants and components, only use surgical instruments manufactured by Plenum, so that there is no inconsistency between the fittings and irreversible damage to the implants and components.
- The list of surgical instruments can be found in the respective printed or electronic communication materials or on the Plenum website (www.plenum.bio/ifu). The use of Plenum instruments must follow the manufacturer's instructions;
- Surgical and prosthetic instruments are not the main object of these instructions for use, i.e. they have separate registrations with ANVISA. They must therefore be purchased separately.

14. STERILIZATION

- The Healing abutment Ti is supplied in sterile condition, sterilized by Gamma Irradiation;
- Product for single use only and cannot be reprocessed, resterilized or reused;
- Before use, check the integrity of the packaging and the expiration date. Do not use the product if the packaging is broken or has expired;
- Date of manufacture, expiry date, code, batch number and other information regarding identification and traceability are on the product label;
- To maintain sterility, the product should only be opened at the time of use and handled in appropriate, aseptic environments;
- Opening the packaging to use the product should only be done by qualified and trained personnel.

15. SHELF LIFE OF THE PRODUCT

- Expiration date indicated on the label.
- Before use, check the expiration date. Do not use the product if it has expired.

16. IDENTIFICATION AND TRACEABILITY

In order to guarantee the traceability of the Healing abutment Ti and comply with the requirements of current legislation, the surgeon or someone on their team must keep information on the membrane used in the patient's medical record. This information must also be passed on to the product distributor and the patient, in order to complete the traceability cycle.

It is the responsibility of the surgeon in charge and/or the hospital/healthcare facility to ensure complete traceability of the Ti Scar, by writing down the code and batch of the product used on the patient's medical record, or by identifying it using one of the 3 traceability labels that come with the packaging, and attaching this information to the patient's medical record in order to maintain traceability. In addition, this information should be provided to the patient so that they have information about the product used in their procedure.

The labels contain the following information necessary for the traceability of the product used:

- Manufacturer identification;
- Name or commercial model of the product;
- Product code;
- Product batch number;
- ANVISA registration number.

Traceability information is necessary for notification by the health service and/or the patient themselves to the National Health Surveillance Agency (ANVISA) and the manufacturer, when adverse events occur, in order to conduct the appropriate investigations.

17. CONDITIONS STORAGE STORAGE CONDITIONS, TRANSPORTATION, CONSERVATION AND HANDLING

- The packaging must be intact upon receipt; do not use the product if the packaging is damaged or broken, or if the expiration date has passed. If the packaging is damaged, the product should be considered non-sterile and discarded;
- The product must be stored and transported in a clean environment, protected from sources of heat or humidity, away from sunlight and free from the effects of the weather, so that prevent any damage or alteration to the packaging or its physical form that makes it impossible to use;
- The product should be stored on shelves or in drawers/cabinets in such a way that it is easy for the operator to see and handle. Do not store the product near light bulbs, so as not to affect the information on the packaging;
- The product should only be opened at the time of use in a sterile surgical field;
- The product should only be handled by qualified and trained personnel.

18. DISPOSAL OF THE PRODUCT

- Discard the product if any of the successive packages ensuring sterility are damaged or opened;
- After use or any removal of the product from the patient, the Healing abutment Ti that is not subject to further analysis or study must be disposed of by the health service establishment (hospitals, clinics and practices);
- It is recommended that the product be completely disposed of so that it cannot be reused.

The disposal of the product is the sole responsibility of the institution/health service establishment;

- After de-characterization, the product must be identified as unfit for use and disposed of in accordance with health service waste management standards (hospitals, practices and clinics), in accordance with applicable environmental legislation;
- Healing abutment Ti that is rendered unusable during surgery should also be disposed of following the guidelines described above.

19. CUSTOMER SERVICE

- If you would like to obtain the Instructions for Use in printed form and at no additional cost (including postage), please call (11) 3109-9000, e-mail (relacionamento@plenum.bio) or visit the website (www.plenum.bio/ifu);
- If the product presents any adverse event, with risk potential that could harm the patient, generates or has the potential to harm or threaten public health, is out of its specifications, or is generating any dissatisfaction, the healthcare professional must generate a technical report containing the details of the non-conformity, with all the complainant's data, product code and batch, and patient data; notifying Plenum's Customer Relations Service directly; - If necessary, send the product clean, decontaminated, sterilized, packaged and duly identified, with a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logístico - CEP 13212-213 - Jundiaí - SP - Brazil. Telephone: (11) 3109-9000; E-mail: relacionamento@plenum.bio - ANVISA also provides a channel for receiving notifications of incidents, adverse events and technical complaints related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at www.anvisa.gov.br, link NOTIVISA.

20. LIMITED WARRANTY

- Products manufactured by Plenum have a factory guarantee against manufacturing defects. Defects in the product must be reported immediately to the manufacturer, observing the legal deadline, attaching a technical report to the product detailing the non-conformity found. This report must be generated and signed by a qualified healthcare professional. In the event of a risk to the patient's health, the relevant ANVISA regulations must also be observed (see Customer Service section).

The warranty on products manufactured by Plenum is subject to strict observance and compliance with the indications for use, contraindications, precautions, warnings and other information contained in the instructions for use. Failure to do so will release the manufacturer and distributor of the products from any liability whatsoever, especially in the following cases:

- Reprocessing, resterilization or reuse of the product by the user;
- Inappropriate use or use other than in a surgical environment;
- If the product has not been used in full compliance with the instructions for use.

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