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E-mail: relacionamento@plenum.bio

Site: www.plenum.bio

INSTRUCTIONS FOR USE LABORATORY COMPONENT

INSTRUCTIONS FOR USE

TECHNICAL NAME: Not applicable. Product for laboratory use.

COMMERCIAL NAME: Laboratory component.

MODELS: Detailed product description see table 1.

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Read all Instructions for Use carefully before use. Follow all directions, warnings, and precautions mentioned in these Instructions for Use. Failure to observe these points may lead to complications. Product for Dental Use. Product Supplied Non-Sterile. Indefinite Shelf Life. Manufacturer recommends single use. Do not reuse.

ANVISA Registration No.: Not applicable. Product for laboratory use. Code, Lot No., and Date of Manufacture: See product labeling.

Rev. 03

SUMMARY




1. DETAILED DESCRIPTION OF THE PRODUCT	3
2. PRODUCT COMPOSITION	4
3. DEVICE PRESENTATION FORM	4
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	4
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	4
6. INSTRUCTIONS FOR USE	4
7. PRECAUTIONS, RESTRICTIONS, AND WARNING	5
8. CONTRAINDICATIONS	5
9. POSSIBLE ADVERSE EFFECTS	5
10. DESCRIPTIONS SYMBOLS	6
11. VALIDITY OF THE LABORATORY PRODUCT	6
12. IDENTIFICATION AND TRACEABILITY	6
13. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS	6
14. DISPOSAL OF DEVICE	7
15. CUSTOMER SERVICE	7
16. LIMITED WARRANTY	7

1. DETAILED DESCRIPTION OF THE PRODUCT

The Laboratory Component is a dental laboratory product available in the Hybrid Analog model, as listed in Table 1.

The Laboratory Component – Hybrid Analog Model is specifically designed to be fixed to a resin or plaster model to replicate the positioning of the implant or prosthetic component that has been installed in a patient. Table 1 below shows an illustrative image, codes, models, and main dimensions of the Hybrid Analog.

TABLE 1. Illustrative images of the Hybrid Analog - codes, models, and dimensions.

Illustrative image	Code	Model	Main Dimensions (mm)
	018-21-C26-00-00	HYBRID - PLENUM IMPLANT SL	Ø3.0
	018-22-C26-00-00	HYBRID - PLENUM IMPLANT RE	Ø3.5
	018-23-C26-00-00	HYBRID - PLENUM IMPLANT SH	Ø4.5
	018-20-C26-20-00	HYBRID MULTI ABUTMENT	Ø4.5
	018-20-C26-20-36	MULTI ABUTMENT Ø3.6	Ø3.6
	018-20-C26-90-35	ANALOG UNIV.Ø3.5 H4	Ø3.5 x 4
	018-20-C26-91-35	ANALOG UNIV.Ø3.5 H6	Ø3.5 x 6
	018-20-C26-90-35	ANALOG UNIV.Ø4.5 H4	Ø4.5 x 4
	018-20-C26-91-35	ANALOG UNIV.Ø4.5 H6	Ø4.5 x 6
	018-20-C26-50-33	HYBRID PLENUM ABT. Ø3.5X4MM	Ø3.5 x 4
	018-20-C26-51-33	HYBRID PLENUM ABT. Ø3.5X6MM	Ø3.5 x 6
	018-20-C26-50-45	HYBRID PLENUM ABT. Ø4.5X4MM	Ø4.5 x 4
	018-20-C26-51-45	HYBRID PLENUM ABT. Ø4.5X6MM	Ø4.5 x 6
	018-20-C26-20-39	HYBRID TI-BASE Ø3.9	Ø3.9
	018-20-C26-20-45	HYBRID TI-BASE Ø4.5	Ø4.5
	018-20-C26-15-00	HYBRID EQUATOR ABUTMENT	Ø4.5
	018-25-C26-00-00	HYBRID - PLENUM IMPLANT MF	Ø4
	018-24-C26-00-00	HYBRID - PLENUM IMPLANT UF	Ø3.5

ACCESSORIES:

Illustrative image	Code	Model	Main Dimensions (mm)
	018-15-F20-00-34	HYBRID ANALOG DRIVER	Ø1.4 / Ø1.6

2. PRODUCT COMPOSITION

The Laboratory Component, in the models described in Table 1, is manufactured in Grade 4 Pure Titanium, in accordance with ASTM F67 and ISO 5832-2 specifications.

The accessory is manufactured in Stainless Steel in accordance with ASTM F899 specifications.

3. DEVICE PRESENTATION FORM

The Hybrid Analog has an insertion key (accessory) that must be used to assist in installing the analog in a resin model.

The Laboratory Component is supplied individually packaged, in a non-sterile condition, and presented as follows:

- Primary packaging in the form of polyethylene terephthalate (PET) blister pack and heat-sealed with surgical-grade Tyvek paper;
- Outer packaging in the form of a rigid cardboard box containing the primary packaging and three identification and traceability labels;
- If you would like to receive the Instructions for Use in printed format at no additional cost (including shipping), please request them by phone at (11) 3109-9000, by email at (relacionamento@plenum.bio), or through the website (www.plenum.bio);
- See below for illustrative images of the primary and external packaging of a Laboratory Component model in the form in which it will be sold.

FIGURE 1. Primary Packaging of the Laboratory Component

FIGURE 2. External Packaging of the Laboratory Component.

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

Laboratory Components - Hybrid Analog model are dental products for laboratory use, designed to accurately transfer the positioning of implants or prosthetic components to a resin or plaster model for clinical and/or laboratory use.

5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

Laboratory Components - Hybrid Analog models function as auxiliary elements in the manufacture of dental prostheses to be installed on Plenum Dental Implants.

Through the use of Hybrid Analogs, it is possible to transfer the positioning of the implant or prosthetic component to a working model made by 3D printing or plaster, so that dental prostheses can be adapted and finalized before installation in the mouth.

In the virtual planning, it is possible to select the prosthetic components such as Plenum Abutment, Ti-Base, and Mini-Abutment to be used later on the Hybrid Analog, allowing greater predictability in the prosthetic results.

6. INSTRUCTIONS FOR USE

- Read all instructions for use, recommendations, and guidelines carefully before using the Laboratory Component.
- Before use, check the product packaging. Do not use the component if the packaging is tampered with or damaged;
- Product supplied in non-sterile condition.
- The Laboratory Component family consists of Hybrid Analog models, whose main function is to be used as auxiliary elements in the manufacture of dental prostheses on implants. (see specific functions in the topic “Detailed Product Description”), and it is the responsibility of the dental surgeon to choose the product models with the most appropriate characteristics and dimensions, after individual clinical analysis of each clinical case;
- The techniques for using and applying the component models that are part of this family vary according to the preference of the dental surgeon and prosthetic technician, who is responsible for choosing the best clinical approach;
- Carefully observe the correct fit between the Laboratory Component and the working model, whether made of plaster or resin;
- For proper installation of the prosthetic component, it is essential to use Plenum instruments, according to the protocol recommended by Plenum (see information on the website: www.plenum.bio);
- Before starting the procedure, make sure that the necessary components, models, and instruments are intact and complete;
- The components to be used must be chosen by the surgeon, dentist, who should consider, among other things, their previous clinical experience and the techniques established in scientific literature;
- Care must also be taken with the contact surface between the component and the working model. Both surfaces must be free of particles and solid materials that could prevent the product from fitting and fixing perfectly;

- In order to protect the patient and their surgeon, it is important to ensure the correct and complete traceability of the component used by noting the product code and batch number in the patient's protocol.

Handling and Manipulation of the Laboratory Component:

1) The dental surgeon or prosthodontist must select the Laboratory Component – Hybrid Analog Model according to the implant installed or intermediate component, if used.

2) After selecting the Laboratory Component, installation should be performed with adaptation to a plaster or resin model.

The fit must be checked carefully by the dentist and/or prosthodontist.

NOTE: For use with the resin-printed model, the Analog should be secured with the key, rotating the assembly until the Analog indexes on the printed model, thus preventing rotation. It should be pressed down until you feel the Analog lock or click into place.

The Analog key must be unscrewed from the Analog before removal.

3) The Laboratory Components are available in a virtual library that is pre-installed in the planning software. At this stage, the Plenum Abutments, Ti-Bases, and Mini-Abutments can be selected according to the needs of the clinical case.

4) Once planning is complete, the digital files for the prosthesis and model to be manufactured will be obtained for sending to production.

5) The prostheses will be manufactured by machining using specific milling machines according to the selected material.

Resin models can be manufactured by additive manufacturing using 3D printers so that the prostheses can be finished before installation in the patient.

7. PRECAUTIONS, RESTRICTIONS, AND WARNING

- The dental surgeon and/or prosthetist should not use the Laboratory Component before carefully reading all instructions for use and information contained in the packaging. All information on use, warnings, precautions, and recommendations mentioned in these instructions for use must be observed. Failure to observe these points may lead to complications for the dental surgeon and, above all, for the patient.

- Single-use product, supplied in a non-sterile condition.

- Before use, check the integrity of the packaging. Do not use the product if the packaging is damaged.

- The date of manufacture, product code, batch number, and other information regarding identification and traceability are on the product label.

- The package should only be opened by qualified and trained personnel for product use;

- Before the procedure, the dentist must carry out rigorous clinical planning, selecting the model and size of the prosthetic component most appropriate for the patient's condition;

- For safety and effectiveness and to not compromise the intended use, only prosthetic components and implants manufactured by Plenum should be used. Components and/or implants manufactured by other companies may

cause dimensional incompatibility, improper fit, micro-movements, increase the risk of corrosion (electrolytic or galvanic) and fracture of these products, compromising the outcome of the surgical and prosthetic procedure;

- Before use, special attention should be paid to the insertion keys. Damage to the tip of the keys may prevent the component from being installed correctly, compromising the intended performance. Therefore, do not use this product if any specific tools are damaged.

- The selection and incorrect choice of the component to be used, as well as errors in indication, handling, and installation, may cause damage to the structure of the intermediate component, which may contribute to failure due to fatigue, fracture, or loosening of the component;

- The selection and incorrect choice of the component to be used, as well as errors in indication, handling, and installation, may interfere with virtual planning, resulting in a poorly fitted prosthesis;

- Do not force the introduction of a prosthetic component beyond the established depth. This action may cause deformation of the model, the hexagonal connection of the component and/or the insertion key, and fracture of the product.

- The success of a dental implant and prosthetic component placement procedure is related to the correct indication and selection of products, correct surgical technique, as well as the time and quality of osseointegration, which can be affected by numerous biological, clinical, and biomechanical factors that are extrinsic to product quality. These factors can cause early failure and loss of products and compromise the surgical procedure. Therefore, strict compliance with all indications for use, contraindications, precautions, warnings, and other observations contained in these instructions for use is essential for the success of the procedure.

- The patient must be informed of all limitations, risks, complications, warnings, and possible adverse effects arising from the implant dentistry procedure, with the use of implants and prosthetic components.

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9. POSSIBLE ADVERSE EFFECTS

As this is a product for laboratory use, it has no adverse effects on patients directly related to the use of the Laboratory Components. However, the inappropriate choice of model or size of these components for implantation or use, or the use of inappropriate techniques, may render the procedure for manufacturing a suitable dental prosthesis unfeasible.

10. DESCRIPTIONS SYMBOLS

Table 2 describes the symbols that may be printed on the product labeling (primary and secondary packaging). Please check the physical labeling to see which symbols apply to the product.

TABLE 2. Symbols representing standards and compliance associated with the laboratory component and its use.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Non-sterile
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
 plenum.bio/ifu	Consult electronic instructions for use (plenum.bio/ifu)
Rev xx	IU Revision Number
	Consult the instructions for use in electronic format (QR code)
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
Reg. ANVISA	ANVISA Registration Number
M3 Health Ind. Com. de Prod. Med. Odont. e correlatos S.A. CNPJ:12.568.799/0001-04	Manufacturer's business name

NOTE. The graphic symbols related to the labeling of the Laboratory Component, referenced above, meet the requirements established in ISO 15223-1.

11. VALIDITY OF THE LABORATORY PRODUCT

The Laboratory Component is supplied in a non-sterile condition. This product has an indefinite shelf life in the form in which it was originally marketed.

12. IDENTIFICATION AND TRACEABILITY

The label contains the following information necessary for identifying the component used:

- Manufacturer identification;
- Contents;
- Trade name;
- Commercial model;
- Raw material;
- Date of manufacture;
- Expiration date;
- Product code;
- Product batch;
- Bar code.
- Symbols used on product labeling, as shown below.
- Symbol indicating that the product is not sterile.
- ANVISA registration number: Not applicable. Product for laboratory use.

13. STORAGE, TRANSPORT, CONSERVATION, AND HANDLING CONDITIONS

- The product packaging must be intact upon receipt. Do not use the product if the packaging is damaged or tampered with. If the packaging is damaged, the product must be discarded.
- The identification data is on the product label;
- The product must be stored and transported in a clean environment, protected from sources of heat or moisture, away from sunlight and free from the effects of weather to prevent any damage or alteration

to the packaging or its physical form that would make it impossible to use.

- The product should be stored on shelves or in drawers/cabinets in such a way as to facilitate viewing and handling by the operator. Do not store the product near light bulbs, so that the information on the packaging is not affected.
- The product should only be handled by personnel who are qualified and trained to perform this procedure.

14. DISPOSAL OF DEVICE

- After use, any component that is not subject to further analysis or study must be disposed of by the hospital/health care facility.
- It is recommended that the product be completely decharacterized so that it cannot be reused. The decharacterization of the product is the sole responsibility of the institution/health service establishment;
- After decharacterization, the product must be identified as unsuitable for use and disposed of in accordance with health service waste management standards (hospitals, offices, and clinics) in accordance with applicable environmental legislation;
- Products that are rendered unusable during laboratory processing must also be disposed of in accordance with the guidelines described above.

15. CUSTOMER SERVICE

- If the product presents any adverse event with potential risk that could harm the patient, generate or have the potential to cause injury or threaten public health, is outside its specifications, or is causing any dissatisfaction, the healthcare professional must generate a technical report containing the details of the non-compliance, with all the complainant's data, product code and batch, and patient data; notifying Plenum's Customer Service directly;
- If necessary, send the product clean, decontaminated, packaged, properly identified, and with a description of the nonconformities to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, No. 640 - Lot 15 - Block B - Jardim Ermida I, Multivias II - Industrial and Logistics Hub - CEP 13212-213 - Jundiaí - SP - Brazil. Phone/Fax: (11) 3109-9000; Email: relacionamento@plenum.bio
- ANVISA also provides a channel for receiving notifications of incidents, adverse events, and technical complaints related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency - ANVISA, at the electronic address www.anvisa.gov.br, link NOTIVISA.

16. LIMITED WARRANTY

- Products manufactured by Plenum are covered by a factory warranty against manufacturing defects. Any product defects must be reported immediately to the manufacturer, within the legal deadline, attaching a technical report with details of the non-conformity found. This report must be generated and signed by a qualified and trained healthcare professional, observing the relevant ANVISA standards (see Customer Service above) in case of risk to the patient's health.

The warranty for products manufactured by Plenum is conditional upon strict compliance with and adherence to the instructions for use, contraindications, precautions, warnings, and other information contained in the instructions for use. Failure to comply with and adhere to these instructions exempts the manufacturer and distributor of the products from any liability, particularly in the following cases:

- Reuse of the product by the user;
- Improper use or use other than in a surgical environment;
- If the product has not been used in full compliance with the instructions for use.

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