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PROSTHETIC
COMPONENT
English

COMPONENTE
PROTÉSICO
Español

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E-mail: relacionamento@plenum.bio

Site: www.plenum.bio

INSTRUCTIONS FOR USE PROSTHETIC COMPONENT

INSTRUCTIONS FOR USE

TECHNICAL NAME: Dental Implant Components

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Carefully read these Instructions for Use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for Use. A failure to observe these items may result in complications. Device for dental use. Device supplied without sterilization. Sterilize before use. Expiration date: Indeterminate. Single-use only. Prohibited reuse.

ANVISA DEVICE REGISTRATION N°.: 81684340013

Device Code, batch code, and date of manufacture: See device labeling.

Rev.08

SUMMARY

1. DETAILED PRODUCT DESCRIPTION	3
2. DEVICE COMPOSITION	6
3. DEVICE PRESENTATION FORM	6
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	7
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	7
6. INSTRUCTIONS FOR USE	7
7. PRECAUTIONS	8
8. WARNINGS	9
9. CONTRAINDICATIONS	9
10. POSSIBLE ADVERSE EFFECTS	9
11. MAGNETIC RESONANCE IMAGING (MRI) - INFORMATION	9
12. SYMBOL DESCRIPTIONS	10
13. AUXILIARY INSTRUMENTS	10
14. STERILIZATION	10
15. EXPIRATION DATE	11
16. IDENTIFICATION AND TRACEABILITY	11
17. STORAGE, TRANSPORT, AND HANDLING CONDITIONS	11
18. DISPOSAL OF DEVICE	11
19. CUSTOMER SERVICES	11
20. LIMITED WARRANTY	12

1. DETAILED PRODUCT DESCRIPTION

Prosthetic Components are available in models according to Tables 1 and 2.

The UCLA Component has a cylindrical shape in the form of a pillar, a metal base made of Cobalt-Chromium-Molybdenum (CoCr) alloy, and a cylindrical complement made of calcinable polymer (acrylate). The UCLA is used for the fabrication of permanent dental prostheses, in which the waxing up and subsequent over-molding is done in the prosthetic laboratory. The acrylic complement supports the molding resin, which is fused and replaced by a biocompatible metal alloy (see Product Composition and Acceptable Metal Combinations) and joins the supplied metal base together with UCLA. In this way, an individualized prosthetic component is formed. It has an anti-

rotational connection at its end and is available on the Slim (SL), Regular (RE), Unifit (UF), Morsefit (MF), and Short (SH) platforms, as detailed in Table 1. After over-molding, the UCLA becomes an intermediate prosthetic component, which is mounted and fixed directly on the dental implant by means of a screw on which, in a subsequent procedure, a dental prosthesis will be fixed. Therefore, UCLA is an element that makes the connection between the Plenum Dental Implant installed next to the bone tissue and a dental prosthesis that will be used.

Table 1 presents the illustrative images, codes, models, prosthetic platforms, and the main dimensions of the UCLA Components.

TABLE 1: Illustrative images of UCLA Components, codes, models, and dimensions.

Imagem ilustrativa	Code*	Trade Name / Model	Composition of the trade model *	Diameter (mm)	Length (mm)	Prosthetic Interface
	018-21-C10-33-00	UCLA SL Ø3.5X0MM	UCLA SL Ø3.5 X 0MM + HEXAGONAL SCREW SL 0MM	3.5	0.0	SL
	018-21-C10-33-01	UCLA SL Ø3.5X1MM	UCLA SL Ø3.5 X 1MM + HEXAGONAL SCREW SL 1MM	3.5	1.0	SL
	018-21-C10-33-02	UCLA SL Ø3.5X2MM	UCLA SL Ø3.5 X 2MM + HEXAGONAL SCREW SL 2MM	3.5	2.0	SL
	018-21-C10-33-03	UCLA SL Ø3.5X3MM	UCLA SL Ø3.5 X 3MM + HEXAGONAL SCREW SL 3MM	3.5	3.0	SL
	018-21-C10-33-04	UCLA SL Ø3.5X4MM	UCLA SL Ø3.5 X 4MM + HEXAGONAL SCREW SL 4MM	3.5	4.0	SL
	018-21-C10-33-05	UCLA SL Ø3.5X5MM	UCLA SL Ø3.5 X 5MM + HEXAGONAL SCREW SL 5MM	3.5	5.0	SL
	018-21-C10-33-06	UCLA SL Ø3.5X6MM	UCLA SL Ø3.5 X 6MM + HEXAGONAL SCREW SL 6MM	3.5	6.0	SL
	018-22-C10-33-00	UCLA RE Ø3.6X0MM	UCLA RE Ø3.6 X 0MM + HEXAGONAL SCREW RE 0MM	3.6	0.0	RE
	018-22-C10-33-01	UCLA RE Ø3.6X1MM	UCLA RE Ø3.6 X 1MM + HEXAGONAL SCREW RE 1MM	3.6	1.0	RE
	018-22-C10-33-02	UCLA RE Ø3.6X2MM	UCLA RE Ø3.6 X 2MM + HEXAGONAL SCREW RE 2MM	3.6	2.0	RE
	018-22-C10-33-03	UCLA RE Ø3.6X3MM	UCLA RE Ø3.6 X 3MM + HEXAGONAL SCREW RE 3MM	3.6	3.0	RE
	018-22-C10-33-04	UCLA RE Ø3.6X4MM	UCLA RE Ø3.6 X 4MM + HEXAGONAL SCREW RE 4MM	3.6	4.0	RE
	018-22-C10-33-05	UCLA RE Ø3.6X5MM	UCLA RE Ø3.6 X 5MM + HEXAGONAL SCREW RE 5MM	3.6	5.0	RE
	018-22-C10-33-06	UCLA RE Ø3.6X6MM	UCLA RE Ø3.6 X 6MM + HEXAGONAL SCREW RE 6MM	3.6	6.0	RE
	018-22-C10-45-00	UCLA RE Ø4.5X0MM	UCLA RE Ø4.5 X 0MM + HEXAGONAL SCREW RE 0MM	4.5	1.0	RE
	018-22-C10-45-01	UCLA RE Ø4.5X1MM	UCLA RE Ø4.5 X 1MM + HEXAGONAL SCREW RE 1MM	4.5	2.0	RE
	018-22-C10-45-02	UCLA RE Ø4.5X2MM	UCLA RE Ø4.5 X 2MM + HEXAGONAL SCREW RE 2MM	4.5	3.0	RE
	018-22-C10-45-03	UCLA RE Ø4.5X3MM	UCLA RE Ø4.5 X 3MM + HEXAGONAL SCREW RE 3MM	4.5	4.0	RE
	018-22-C10-45-04	UCLA RE Ø4.5X4MM	UCLA RE Ø4.5 X 4MM + HEXAGONAL SCREW RE 4MM	4.5	5.0	RE
	018-22-C10-45-05	UCLA RE Ø4.5X5MM	UCLA RE Ø4.5 X 5MM + HEXAGONAL SCREW RE 5MM	4.5	6.0	RE
	018-22-C10-45-06	UCLA RE Ø4.5X6MM	UCLA RE Ø4.5 X 6MM + HEXAGONAL SCREW RE 6MM	4.5	1.0	RE
	018-23-C10-33-00	UCLA SH Ø4.5X0MM	UCLA SH Ø4.5 X 0MM + HEXAGONAL SCREW SH 0MM	4.5	0.0	SH
	018-23-C10-33-01	UCLA SH Ø4.5X1MM	UCLA SH Ø4.5 X 1MM + HEXAGONAL SCREW SH 1MM	4.5	1.0	SH
	018-23-C10-33-02	UCLA SH Ø4.5X2MM	UCLA SH Ø4.5 X 2MM + HEXAGONAL SCREW SH 2MM	4.5	2.0	SH
	018-23-C10-33-03	UCLA SH Ø4.5X3MM	UCLA SH Ø4.5 X 3MM + HEXAGONAL SCREW SH 3MM	4.5	3.0	SH
	018-23-C10-33-04	UCLA SH Ø4.5X4MM	UCLA SH Ø4.5 X 4MM + HEXAGONAL SCREW SH 4MM	4.5	4.0	SH
	018-23-C10-33-05	UCLA SH Ø4.5X5MM	UCLA SH Ø4.5 X 5MM + HEXAGONAL SCREW SH 5MM	4.5	5.0	SH
	018-23-C10-33-06	UCLA SH Ø4.5X6MM	UCLA SH Ø4.5 X 6MM + HEXAGONAL SCREW SH 6MM	4.5	6.0	SH
	018-24-C10-33-01	UCLA UF Ø3,6X1MM	UCLA UF Ø3,6 X 1MM + HEXAGONAL SCREW UF 1MM	3.6	1.0	UF
	018-24-C10-33-02	UCLA UF Ø3,6X2MM	UCLA UF Ø3,6 X 2MM + HEXAGONAL SCREW UF 2MM	3.6	2.0	UF
	018-24-C10-33-03	UCLA UF Ø3,6X3MM	UCLA UF Ø3,6 X 3MM + HEXAGONAL SCREW UF 3MM	3.6	3.0	UF

	018-24-C10-33-04	UCLA UF Ø3,6X4MM	UCLA UF Ø3,6 X 4MM + HEXAGONAL SCREW UF 4MM	3.6	4.0	UF
	018-24-C10-33-05	UCLA UF Ø3,6X5MM	UCLA UF Ø3,6 X 5MM + HEXAGONAL SCREW UF 5MM	3.6	5.0	UF
	018-24-C10-33-06	UCLA UF Ø3,6X6MM	UCLA UF Ø3,6 X 6MM + HEXAGONAL SCREW UF 6MM	3.6	6.0	UF
	018-24-C10-45-01	UCLA UF Ø4,5X1MM	UCLA UF Ø4,5 X 1MM + HEXAGONAL SCREW UF 1MM	4.5	1.0	UF
	018-24-C10-45-02	UCLA UF Ø4,5X2MM	UCLA UF Ø4,5 X 2MM + HEXAGONAL SCREW UF 2MM	4.5	2.0	UF
	018-24-C10-45-03	UCLA UF Ø4,5X3MM	UCLA UF Ø4,5 X 3MM + HEXAGONAL SCREW UF 3MM	4.5	3.0	UF
	018-24-C10-45-04	UCLA UF Ø4,5X4MM	UCLA UF Ø4,5 X 4MM + HEXAGONAL SCREW UF 4MM	4.5	4.0	UF
	018-24-C10-45-05	UCLA UF Ø4,5X5MM	UCLA UF Ø4,5 X 5MM + HEXAGONAL SCREW UF 5MM	4.5	5.0	UF
	018-24-C10-45-06	UCLA UF Ø4,5X6MM	UCLA UF Ø4,5 X 6MM + HEXAGONAL SCREW UF 6MM	4.5	6.0	UF
	018-25-C10-33-01	UCLA MF Ø3,6 X 1MM	UCLA MF Ø3,6 X 1MM + HEXAGONAL SCREW MF 1MM	3.6	1.0	MF
	018-25-C10-33-02	UCLA MF Ø3,6 X 2MM	UCLA MF Ø3,6 X 2MM + HEXAGONAL SCREW MF 2MM	3.6	2.0	MF
	018-25-C10-33-03	UCLA MF Ø3,6 X 3MM	UCLA MF Ø3,6 X 3MM + HEXAGONAL SCREW MF 3MM	3.6	3.0	MF
	018-25-C10-33-04	UCLA MF Ø3,6 X 4MM	UCLA MF Ø3,6 X 4MM + HEXAGONAL SCREW MF 4MM	3.6	4.0	MF
	018-25-C10-33-05	UCLA MF Ø3,6 X 5MM	UCLA MF Ø3,6 X 5MM + HEXAGONAL SCREW MF 5MM	3.6	5.0	MF
	018-25-C10-33-06	UCLA MF Ø3,6 X 6MM	UCLA MF Ø3,6 X 6MM + HEXAGONAL SCREW MF 6MM	3.6	6.0	MF
	018-25-C10-45-01	UCLA MF Ø4,5 X 1MM	UCLA MF Ø4,5 X 1MM + HEXAGONAL SCREW MF 1MM	4.5	1.0	MF
	018-25-C10-45-02	UCLA MF Ø4,5 X 2MM	UCLA MF Ø4,5 X 2MM + HEXAGONAL SCREW MF 2MM	4.5	2.0	MF
	018-25-C10-45-03	UCLA MF Ø4,5 X 3MM	UCLA MF Ø4,5 X 3MM + HEXAGONAL SCREW MF 3MM	4.5	3.0	MF
	018-25-C10-45-04	UCLA MF Ø4,5 X 4MM	UCLA MF Ø4,5 X 4MM + HEXAGONAL SCREW MF 4MM	4.5	4.0	MF
	018-25-C10-45-05	UCLA MF Ø4,5 X 5MM	UCLA MF Ø4,5 X 5MM + HEXAGONAL SCREW MF 5MM	4.5	5.0	MF
	018-25-C10-45-06	UCLA MF Ø4,5 X 6MM	UCLA MF Ø4,5 X 6MM + HEXAGONAL SCREW MF 6MM	4.5	6.0	MF
	018-21-C11-33-00	UCLA SMART SL Ø3.5X0MM	UCLA SMART SL Ø3.5 X 0MM + HEXALOBULAR SCREW SL 0MM	3.5	0.0	SL
	018-21-C11-33-01	UCLA SMART SL Ø3.5X1MM	UCLA SMART SL Ø3.5 X 1MM + HEXALOBULAR SCREW SL 1MM	3.5	1.0	SL
	018-21-C11-33-02	UCLA SMART SL Ø3.5X2MM	UCLA SMART SL Ø3.5 X 2MM + HEXALOBULAR SCREW SL 2MM	3.5	2.0	SL
	018-21-C11-33-03	UCLA SMART SL Ø3.5X3MM	UCLA SMART SL Ø3.5 X 3MM + HEXALOBULAR SCREW SL 3MM	3.5	3.0	SL
	018-21-C11-33-04	UCLA SMART SL Ø3.5X4MM	UCLA SMART SL Ø3.5 X 4MM + HEXALOBULAR SCREW SL 4MM	3.5	4.0	SL
	018-21-C11-33-05	UCLA SMART SL Ø3.5X5MM	UCLA SMART SL Ø3.5 X 5MM + HEXALOBULAR SCREW SL 5MM	3.5	5.0	SL
	018-21-C11-33-06	UCLA SMART SL Ø3.5X6MM	UCLA SMART SL Ø3.5 X 6MM + HEXALOBULAR SCREW SL 6MM	3.5	6.0	SL
	018-22-C11-33-00	UCLA SMART RE Ø3.6X0MM	UCLA SMART RE Ø3.6 X 0MM + HEXALOBULAR SCREW RE 0MM	3.6	0.0	RE
	018-22-C11-33-01	UCLA SMART RE Ø3.6X1MM	UCLA SMART RE Ø3.6 X 1MM + HEXALOBULAR SCREW RE 1MM	3.6	1.0	RE
	018-22-C11-33-02	UCLA SMART RE Ø3.6X2MM	UCLA SMART RE Ø3.6 X 2MM + HEXALOBULAR SCREW RE 2MM	3.6	2.0	RE
	018-22-C11-33-03	UCLA SMART RE Ø3.6X3MM	UCLA SMART RE Ø3.6 X 3MM + HEXALOBULAR SCREW RE 3MM	3.6	3.0	RE
	018-22-C11-33-04	UCLA SMART RE Ø3.6X4MM	UCLA SMART RE Ø3.6 X 4MM + HEXALOBULAR SCREW RE 4MM	3.6	4.0	RE
	018-22-C11-33-05	UCLA SMART RE Ø3.6X5MM	UCLA SMART RE Ø3.6 X 5MM + HEXALOBULAR SCREW RE 5MM	3.6	5.0	RE
	018-22-C11-33-06	UCLA SMART RE Ø3.6X6MM	UCLA SMART RE Ø3.6 X 6MM + HEXALOBULAR SCREW RE 6MM	3.6	6.0	RE
	018-22-C11-45-00	UCLA SMART RE Ø4.5X0MM	UCLA SMART RE Ø4.5 X 0MM + HEXALOBULAR SCREW RE 0MM	4.5	0.0	RE
	018-22-C11-45-01	UCLA SMART RE Ø4.5X1MM	UCLA SMART RE Ø4.5 X 1MM + HEXALOBULAR SCREW RE 1MM	4.5	1.0	RE
	018-22-C11-45-02	UCLA SMART RE Ø4.5X2MM	UCLA SMART RE Ø4.5 X 2MM + HEXALOBULAR SCREW RE 2MM	4.5	2.0	RE
	018-22-C11-45-03	UCLA SMART RE Ø4.5X3MM	UCLA SMART RE Ø4.5 X 3MM + HEXALOBULAR SCREW RE 3MM	4.5	3.0	RE
	018-22-C11-45-04	UCLA SMART RE Ø4.5X4MM	UCLA SMART RE Ø4.5 X 4MM + HEXALOBULAR SCREW RE 4MM	4.5	4.0	RE
	018-22-C11-45-05	UCLA SMART RE Ø4.5X5MM	UCLA SMART RE Ø4.5 X 5MM + HEXALOBULAR SCREW RE 5MM	4.5	5.0	RE
	018-22-C11-45-06	UCLA SMART RE Ø4.5X6MM	UCLA SMART RE Ø4.5 X 6MM + HEXALOBULAR SCREW RE 6MM	4.5	6.0	RE
	018-23-C11-33-00	UCLA SMART SH Ø4.5X0MM	UCLA SMART SH Ø4.5 X 0MM + HEXALOBULAR SCREW SH 0MM	4.5	0.0	SH
	018-23-C11-33-01	UCLA SMART SH Ø4.5X1MM	UCLA SMART SH Ø4.5 X 1MM + HEXALOBULAR SCREW SH 1MM	4.5	1.0	SH
	018-23-C11-33-02	UCLA SMART SH Ø4.5X2MM	UCLA SMART SH Ø4.5 X 2MM + HEXALOBULAR SCREW SH 2MM	4.5	2.0	SH
	018-23-C11-33-03	UCLA SMART SH Ø4.5X3MM	UCLA SMART SH Ø4.5 X 3MM + HEXALOBULAR SCREW SH 3MM	4.5	3.0	SH
	018-23-C11-33-04	UCLA SMART SH Ø4.5X4MM	UCLA SMART SH Ø4.5 X 4MM + HEXALOBULAR SCREW SH 4MM	4.5	4.0	SH
	018-23-C11-33-05	UCLA SMART SH Ø4.5X5MM	UCLA SMART SH Ø4.5 X 5MM + HEXALOBULAR SCREW SH 5MM	4.5	5.0	SH
	018-23-C11-33-06	UCLA SMART SH Ø4.5X6MM	UCLA SMART SH Ø4.5 X 6MM + HEXALOBULAR SCREW SH 6MM	4.5	6.0	SH


IMPORTANT: * Each commercial product code is composed of the UCLA component model and its respective screw (according to the commercial model composition column).UCLA prosthetic components are compatible with all Plenum Implants. provided that the prosthetic interface corresponds to the internal connection of the implant (SL, RE,SH, UF and MF). UCLA SMART – UCLA model with a side opening of the fixation driver. Length (mm): refers to the height of the transmucosal represented by the letter A in the UCLA technical drawing. For the total length of the product (UCLA + complement), 10.8mm is added the height variables of the letter A.]

The Coping Component has a cylindrical shape with a metallic base of CoCr alloy and a complement of calcinable polymer (acrylate). Coping is used for the preparation of waxing and subsequent over-molding in a prosthetic laboratory, where the acrylic complement, which serves as a support for the molding resin, is fused and replaced by a biocompatible metal alloy (see Product Composition and Acceptable Metal Combinations); thus, joining the metallic base provided together with Coping forming a prosthetic component. After over-molding, Coping becomes an intermediate prosthetic component, which is mounted and fixed directly on the prosthetic abutment by means of

a screw, on which a dental prosthesis will be fixed. Therefore, Coping is an element that makes the connection between the prosthetic abutment, and the dental prosthesis will be used. One of its models has an anti-rotational connection (AR) at its end, according to the details and dimensions described in Table 2.

NOTE: It is worth pointing out that the screw is an integral and exclusive part of the Prosthetic Component models, and that without this item the product will not perform its intended use.

TABLE 2: Illustrative images of Coping Components, codes, models and dimensions.

Illustrative Image	Code*	Trade Name / Model	Composition of the trade model *	Diameter (mm)	Length (mm)
	018-20-C45-20-01	UNIV. ABT. COPING Ø3,5 H4	3.5	8.6	Not applicable
	018-20-C45-20-02	UNIV. ABT. COPING Ø3,5 H6	3.5	10.6	Not applicable
	018-20-C46-20-01	UNIV. ABT. COPING Ø3,5 H4 AR	3.5	8.6	Not applicable
	018-20-C46-20-02	UNIV. ABT. COPING Ø3,5 H6 AR	3.5	10.6	Not applicable
	018-20-C45-20-03	UNIV. ABT. COPING Ø4,5 H4	4.5	8.6	Not applicable
	018-20-C45-20-04	UNIV. ABT. COPING Ø4,5 H6	4.5	10.6	Not applicable
	018-20-C46-20-03	UNIV. ABT. COPING Ø4,5 H4 AR	4.5	8.6	Not applicable
	018-20-C46-20-04	UNIV. ABT. COPING Ø4,5 H6 AR	4.5	10.6	Not applicable
	018-20-C46-21-36	MULTI ABUTMENT COPING Ø3,6 AR	3.6	10.1	Not applicable
	018-20-C45-20-36	MULTI ABUTMENT COPING Ø3,6	3.6	10.1	Not applicable
	018-20-C46-21-48	MULTI ABUTMENT COPING AR	4.8	10.1	Not applicable
	018-20-C45-20-00	MULTI ABUTMENT COPING	4.8	10.1	Not applicable
	018-20-C45-50-33	ABUTMENT COPING Ø3.5 H4	3.5	8.6	Not applicable
	018-20-C45-51-33	ABUTMENT COPING Ø3.5 H6	3.5	10.6	Not applicable
	018-20-C46-50-33	ABUTMENT COPING Ø3.5 H4-AR	3.5	8.6	Not applicable
	018-20-C46-51-33	ABUTMENT COPING Ø3.5 H6-AR	3.5	10.6	Not applicable
	018-20-C47-50-33	ABUTMENT COPING Ø3.5 H4/SCREW	3.5	8.6	Not applicable
	018-20-C47-51-33	ABUTMENT COPING Ø3.5 H6/SCREW	3.5	10.6	Not applicable
	018-20-C48-50-33	ABUTMENT COPING Ø3.5 H4-AR/S	3.5	8.6	Not applicable
	018-20-C48-51-33	ABUTMENT COPING Ø3.5 H6-AR/S	3.5	10.6	Not applicable
	018-20-C45-50-45	ABUTMENT COPING Ø4.5 H4	4.5	8.6	Not applicable
	018-20-C45-51-45	ABUTMENT COPING Ø4.5 H6	4.5	10.6	Not applicable
	018-20-C46-50-45	ABUTMENT COPING Ø4.5 H4-AR	4.5	8.6	Not applicable
	018-20-C46-51-45	ABUTMENT COPING Ø4.5 H6-AR	4.5	10.6	Not applicable
	018-20-C47-50-45	ABUTMENT COPING Ø4.5 H4/SCREW	4.5	8.6	Not applicable
	018-20-C47-51-45	ABUTMENT COPING Ø4.5 H6/SCREW	4.5	10.6	Not applicable
	018-20-C48-50-45	ABUTMENT COPING Ø4.5 H4-AR/S	4.5	8.6	Not applicable
	018-20-C48-51-45	ABUTMENT COPING Ø4.5 H6-AR/S	4.5	10.6	Not applicable

IMPORTANT:

Prosthetic Components - COPING are compatible with their respective Abutments. Multi Abutment Coping with Multi Abutment Plenum, Coping H4 Abutment with Plenum H4 Abutment and Coping H6 Abutment with Plenum H6 Abutment. (Abutments are not part of this family.

COPING AR - Coping with internal hexagon for single prosthesis, when not AR Coping for multiple prosthesis.

COPING SCREW - Coping fixed by means of screw, when without a screw, the Coping is fixed with cement.

Total length (mm) - refers to the sum of the length of the coping + its complement.

2. DEVICE COMPOSITION

The Prosthetic Components have the following integral parts according to the schematic models described in Figures 1, 2 and 3 below:

SCHEMATIC MODEL OF THE PARTS OF THE PROSTHETIC COMPONENT (WITH UNIQUE FIXATION SCREW).


FIGURE 1. UCLA model


FIGURE 2. Multi abutment coping with fixation screw

NOTE: It is worth pointing out that the fixation screw is an integral and exclusive part of the Prosthetic Component models, and that without this item the product will not perform its intended use.


FIGURE 3: Multi abutment coping without fixation screw.

NOTE: For models that do not have a screw, fixation occurs through cementation.

Prosthetic Components are products supplied packed individually, in a non-sterile condition, and manufactured in CoCr alloy (alloy 1 UNS R31537), according to ASTM F1537 and NBR ISO 5832-12, and its complement in acrylic (acrylate).

CHEMICAL COMPOSITION OF COCR ALLOY - ALLOY 1 UNS R31537 (%)

Element/Symbol	Min	Max
Carbon (C)	--	0.14
Chromium (Cr)	26	30
Molybdenum (Mo)	5.0	7.0
Nickel (Ni)	--	1.0
Iron (Fe)	--	0.75
Silicon (Si)	--	1.0
Magnesium (Mg)	--	1.0
Nitrogen (N)	--	0.25
Cobalt (Co)	Chemical balance	

3. DEVICE PRESENTATION FORM

The Prosthetic Components do not have accessories.

The Prosthetic Components in the models, as described in the tables above, are provided packaged unitarily, in non-sterile condition, and presented as follows:

Primary packaging is a polyethylene terephthalate (PET) blister and thermo-sealed with surgical grade paper, containing 01 unit of the device;

External packaging is a carton paper box containing the primary packaging and 03 identification and traceability labels;

It is the responsibility of the hospital/health care facility to ensure complete traceability of the device through the annotation in the patient's chart of the code and lot of the device used, or else the identification through one of the 3 traceability labels that come with the package;

A printed version of the Instructions for Use (including shipping) can be requested at no additional cost at the telephone number (11) 3109-9000, or through the electronic address (relacionamento@plenum.bio), or through the website (www. plenum.bio/ifu);

ATTENTION: The identification of the revision of the instructions for use is informed in the header of this document and corresponds to the current version of the models described in Tables 1 and 2. The user must observe the correlation of the revision of the instructions for use indicated and made available on the website, www.plenum.bio, with the respective device purchased.

See below the illustrative images of a Prosthetic Component model's primary and external packaging in the form in which it will be commercialized.

ILLUSTRATIVE IMAGE OF THE PACKAGING PRIMARY

ILLUSTRATIVE IMAGE OF THE PACKAGING EXTERNAL

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

The Prosthetic Components in the UCLA and Coping models are considered intermediate components. UCLA has the finality of being mounted and fixed directly on the dental implant. While the Coping is to be mounted directly on the prosthetic abutment. On these components, a dental prosthesis will be fixed, that is, they are intermediate elements that make the connection between the dental implant or prosthetic abutment and a dental prosthesis. These prosthetic components, together with the dental prosthesis, sequentially implanted, can be used in single or two-stage processes, single or multiple restorations, cemented or screwed.

5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

Prosthetic Components are dental devices used to wax and over-mold screwed or cemented prosthesis infrastructures installed on intermediates (prosthetic abutment) or directly on the dental implant. Used in the laboratory phase to make the prosthesis on plaster models and analogs of the dental implants. After casting, it becomes part of the final prosthesis to restore the patient's masticatory function and aesthetics in cases of total, partial, and unitary restorations. These prosthetic components function as intermediate elements between the abutment or implant and the dental prosthesis. The fitting and fixation between these prosthetic components and the abutment or implant occur through a screw thread and connection of the Morse Taper type.

6. INSTRUCTIONS FOR USE

Carefully read all instructions for use, recommendations, and guidelines before using the Prosthetic Components;

The family of Prosthetic Components consists of the UCLA and Coping models. The UCLA (see Table 1) is intended to be mounted and fastened directly on the Plenum Dental Implant. While the Coping (see Table 2) is to be mounted directly on the prosthetic abutment;

The techniques of use and application of these prosthetic components vary according to the surgeon's preference, and it is up to him to choose the surgical approach, models and dimensions of the devices, the application technique, as well as the criteria for monitoring and evaluating the results of the surgical and prosthetic procedure;

It is recommended to request diagnostic imaging tests for anatomical visualization of the results obtained (for example periapical radiography and computed tomography);

CLINICAL PHASE:

The Prosthetic Components are restricted to dental use and should only be used by previously trained and qualified dentist, with experience and knowledge of implant dentistry techniques;

Device is supplied in non-sterile condition and must be sterilized before use in case the dental surgeon uses the device in an invasive procedure of transient contact. For more details, see the section "Sterilization";

Before starting clinical and/or surgical procedures, the surgeon must subject

the patient to a rigorous oral antisepsis to prevent the component from coming into contact with contaminating substances;

Before rehabilitating the patient, the dental surgeon must carry out rigorous planning through the use of clinical examinations and images (radiographic and tomographic) to define whether the implants installed are in a position to receive prosthetic components;

Before using the Prosthetic Components, the dentist must strictly verify that the device has the same prosthetic interface as the implant/prosthetic component that is installed;

Before starting the clinical procedure, the dental surgeon must make sure that the prosthetic and auxiliary instructional components are intact and complete;

The surgical procedure, the model and size of the component, the technique, and the installation sequence of the Prosthetic Components must respect the established processes of implant dentistry;

For the use of Prosthetic Components in two-stage procedures, we recommend the previous conditioning of the peri-implant tissues with the use of the Healing Abutment or Cap (this item is not part of this instruction for use, it has a separate ANVISA registration and must be purchased separately);

ATTENTION: The figures are merely illustrative. They do not represent the real dimensions and characteristics of the device.

For the use of Prosthetic Components, use only instruments manufactured by Plenum, which are considered auxiliary instruments, are not part of these instructions for use, and must be purchased separately (see Auxiliary Instruments);

The handling of the device must be carried out carefully, especially during fixing the component on the dental implant. The component must passively sit on the implant platform. Forced fitting, use of excessive force, or oblique insertion must be avoided;

It is highly recommended to use imaging tests to certify the correct placement of the prosthetic component on the implant platform;

Special care must be taken with the hexagonal connection existing on some components and with the driver used for their insertion, as an improper fit between the two elements may cause damage to this connection, compromising device integrity and proper fixation to the implant;

Care must be taken with the surface of the screw thread and the Morse Taper connection in the prosthetic component. Both surfaces must be free of particles and solid materials that can prevent the perfect fitting and fixation next to the implant and dental prosthesis, preventing the proper fixation of these elements;

The dentist must carefully evaluate the introduction and removal of a prosthetic component repeatedly in the same patient, as this can cause deterioration of the screw threads, the Morse Taper, and the connection of the component;

Check the passivity and perform the occlusal and interproximal adjustment after the installation of the prosthesis, avoiding the compromise of the implant/prosthesis set.

LABORATORY PHASE

The Prosthetic Components must be subjected to an over-molding process, where the acrylic complement, which serves as a support for the molding resin, is replaced by a metal alloy.

The laboratory techniques for choosing the most suitable component for over-molding, the choice of alloys, their execution, and dimensioning are the responsibility of the laboratory technician.

UCLA SMART RE

MULTI ABUTMENT COPING

ACCEPTABLE METAL COMPOUNDS:

The relationship of “Acceptable Metal Compounds” was based on ISO 21534. Therefore, there is total compatibility between the raw materials used, and they are in compliance with the manufacture of Prosthetic Components, in addition to the CoCr-based alloys that can be used for over-molding in a prosthetic laboratory.

7. PRECAUTIONS

The dentist should not use this device before carefully reading all instructions and the information in the packaging. Failure to observe these points may lead to complications for the surgery and especially for the patient;

The UCLAS of the SL and RE prosthetic interfaces (Ø 3.5 and 3.6 mm) are indicated for the rehabilitation of the anterior region (canines and incisors);

Date of manufacture, device code, batch number, and other information regarding identification and traceability are on the device labeling;

The opening of the packaging for the use of the device must only be done by qualified and trained personnel to perform this procedure;

Immediate and late postoperative care and the correct oral hygiene technique are extremely important. The patient must be instructed about the risks existing for the healthy maintenance of the implants and components;

The patient should be warned that failure to comply with postoperative instructions (immediate and delayed) may lead to the failure and early loss of the component and implant, requiring new surgery for revision or removal;

For safety and effectiveness and not to compromise the performance of the intended use, only prosthetic components and implants manufactured by Plenum should be used. Components and/ or implants manufactured by other companies can cause dimensional incompatibility, improper fitting, micro-movements, increase the risk of corrosion (electrolytic or galvanic) and fracture of these devices, compromising the result of the surgical and prosthetic procedure;

Before use, special attention should be paid to the insertion drivers. Damage to the tip of the drivers can prevent the correct installation of the component, compromising the intended pro-ethical performance. Therefore, do not use this device if there is any damaged instrument for specific use;

The selection and incorrect choice of the component to be used, as well as errors in the indication, manipulation, and installation, can cause damage to the device’s structure, which can contribute to failure by fatigue, fracture, or loosening;

The use of excessive torque during the insertion of a prosthetic component may lead to device failure. We recommend not exceeding 10 Ncm for fixing the Coping and 20 Ncm for fixing the UCLA;

If a prosthetic component and/or implant is overloaded beyond its functional capacity, excessive bone loss, breakdown of devices and other restoration elements may occur;

Do not force the introduction of a prosthetic component beyond the depth established by the thread boundary. This action may cause deformation in the hexagonal connection of the component and/or the insertion driver and fracture of the device;

Incorrect handling of small sized components within the patient’s mouth may incur a risk of aspiration and/or swallowing of these elements;

The dentist must rigorously evaluate the initial mechanical stability of the prosthetic component. If you have not obtained ideal and safe mechanical stability, the dental surgeon should evaluate and consider a possible removal and/or replacement of the same;

Elderly patients with mental problems, chemically dependent or with any other difficulty in following postoperative recommendations may represent a greater risk of failure of this device, as these patients may ignore the instructions and restrictions arising from this procedure;

Smokers, tobacco users, drug users, alcoholics, diabetics or malnourished patients should be warned about the increased incidence of lack of osseointegration, with consequent failure of implants and components;

The success of a procedure for placing a prosthetic component is related to the correct indication, selection, resistance of the implant and component, correct surgical technique, as well as the quality of osseointegration, which can be affected by numerous biological, clinical, and biomechanical factors, which are external to the quality of the implants and components. These factors can limit the useful life of the devices and consequently cause their failure and early loss. Therefore, faithful compliance with all indications for use, contraindications, precautions, warnings, and other observations contained in these instructions for use are essential for the success of the procedure;

If deformations, scratches, cracks, or dents visible to the component occur, it must be discarded;

The patient must be informed of all limitations, risks, complications, warnings, and possible adverse effects arising from the procedure for installing the dental implant and prosthetic component set;

The patient must be informed of the risks involved if the dental implant and prosthetic component assembly is overloaded beyond its functional capacity, which may cause bone loss, failure of the implant, its components and prosthetic restoration;

The patient must be informed of the risks involved in the application of excessive mechanical load or trauma that affects the operated region, which may imply failure, loss of mechanical stability, loosening, breakage, fracture or loss of the dental implant and its components;

The patient should be adequately oriented regarding the correct oral hygiene and care in the postoperative phase. The patient's ability and willingness to follow the instructions of the dental surgeon is one of the most important aspects of a dental implant placement procedure and its components;

The dental surgeon must clearly guide the patient, who must understand the importance of follow-up in the postoperative period, in which the stability, integrity, and positioning of the prosthetic component and dental implant used must be verified, as well as the osseointegration and biological fixation of the implant in the appropriate time;

In order to protect the patient and his or her surgeon, it is important to ensure the correct and complete traceability of the component used by noting in the patient protocol the device code and batch, or by attaching to the patient protocol the identification and traceability labels that come with the device.

8. WARNINGS

Single-use only, do not resterilize or reuse.

Do not use if package is opened or damaged or if expiration date has expired.

9. CONTRAINDICATIONS

Provisional Components are always used with Plenum Dental Implants (supplied separately). Therefore, the dentist should evaluate the patient's condition prudently before surgery and the use of components and implants, in addition to checking for any contraindication reported below. Contraindications may include, but are not limited to:

An implant that is poorly positioned and/or has insufficient osseointegration to properly support or anchor a dental prosthesis and its components, or has strong restrictions on the use of this device;

Patient with poor bone quality due to osteoporosis or osteopenia, which can cause early failure of the dental implant and component due to lack of osseointegration;

Patients who present signs of allergy or hypersensitivity to the chemical components of the raw materials used in the manufacture of the Prosthetic Components (see "Device Composition");

Insufficient interocclusal or mesio-distal space for prosthetic rehabilitation;

Inflammation and/or periodontal infection;

Recent history of systemic or localized infection;

Parafunctional habits;

Non-treatable occlusion and joint problems;

Use in patients with immunological changes and systemic diseases that result in poor healing or in the absence of integrity of existing tissue at the implant site;

Alcoholic patient or drug user who cannot correctly follow postoperative instructions;

Patient with mental illness or unable to follow the recommendations and care of the postoperative phase;

Patient unable to maintain good oral hygiene;

All contraindications included by the dental-scientific literature known for the use of dental implants and their components.

10. POSSIBLE ADVERSE EFFECTS

No adverse effects are expected as long as the device is used in accordance with its instructions for use.

However, positive results may not be achieved in all surgical and prosthetic procedures. Every procedure presents risks and possibilities of complications, such as infections, allergic drug reactions and anesthetic risks, among others. Possible adverse effects associated with the use of prosthetic components may include, but are not limited to:

Pain or discomfort due to the presence of the prosthetic component and/or dental implant, or due to the surgical procedure;

Lack of osseointegration, reabsorption of bone tissue, or even delay in bone healing, which can result in the failure and breakage of the prosthetic component and/or dental implant;

Loss of primary stability resulting from non-healing of the bone at the appropriate time, or caused by the improper installation of the device;

Fracture, displacement or early or late failure of the component and/or implant due to the technique of improper installation;

Failure, displacement or fracture of the component or dental implant due to excessive loads applied to the devices and/or inadequate bone function;

Loss of the component or implant due to parafunctional habits or poor oral hygiene;

Allergic reaction to the raw material of the component;

Alteration of bone structure due to the imbalance of normal forces/stresses on the bone/implant interface;

Bone necrosis, osteoporosis, insufficient vascularization, bone reabsorption, and insufficient bone formation, all of which can cause the loss of the dental implant fixation and consequently the component;

The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient must be informed by the dental surgeon of all these complications, and of the risks resulting from implant surgery.

11. MAGNETIC RESONANCE IMAGING (MRI) - INFORMATION

The product contains metal parts that may be affected in a magnetic resonance imaging (MRI) environment.

Non-clinical tests and MRI simulations have demonstrated that Plenum products made of metallic material can be used conditionally in MRI (MR Conditional).

A patient with this device can be safely scanned in a magnetic resonance imaging (MRI) system under the following conditions:

Static Magnetic Field: ≤ 3.0 T

Maximum Spatial Gradient Magnetic Field: 4,000 gauss/cm (40 T/m)

Maximum Indication for MRI System - Whole-body human exposure Specific Absorption Rate (SAR): 2.3 W/kg (15 minutes of scanning, normal operating mode);

Under the scanning conditions defined above, the medical device is expected to produce a maximum temperature increase of 0.8 °C after 15 minutes of continuous scanning.

In non-clinical tests, an image artifact originating from the medical device extends approximately 11.3 mm from the device when displayed with a gradient echo pulse sequence in a 3T MRI system.

12. SYMBOL DESCRIPTIONS

Table 3 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

TABLE 3. Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Non-sterile
	Do not resterilize
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Do not re-use
	Consult electronic instructions for use
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)

Proibido reprocessar	Reprocessing prohibited (Brazilian Requirement)
	Conditional use in Magnetic Resonance Imaging

NOTE: The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

13. AUXILIARY INSTRUMENTS

All surgical and prosthetic instruments developed exclusively for the installation of the Prosthetic Component are considered as auxiliary instruments;

Plenum instruments are made of stainless steel, which gives them high strength and durability, according to requirements specified in ASTM-F899. The use of instruments manufactured by Plenum is essential to ensure safety and effectiveness in the installation of components according to the surgical protocol recommended by the manufacturer;

The prosthetic components have a hexagonal-type connection. It is essential that the surgeon use only Plenum drivers to avoid damage to the connection of the components and also to the drivers;

The list of surgical and prosthetic instruments can be found in the respective printed or electronic communication materials or on the Plenum website (www.plenum.bio). The use of Plenum instruments must follow the manufacturer’s instructions;

Surgical and prosthetic instruments are not the main purpose of these instructions for use, that is, they have separate registrations with ANVISA, so they must be purchased separately. Use only instruments manufactured by Plenum so that there is no incongruence between the fittings and irreversible damage to the components and instruments.

14. STERILIZATION

The Prosthetic Components that are part of this family are supplied in non-sterile condition and must be sterilized before use by the hospital/health service establishment;

Prior to sterilization, the devices that are part of this family must be packaged in an appropriate location or in packaging suitable for sterilization, with the aid of clean tongs and/or gloves to avoid contamination. Sterilization packaging is subject to the method of sterilization, transportation and storage, which has a considerable influence on sterilization results. Packaging must be selected in a manner that ensures that devices fit properly inside the packaging. Use a sterilization indicator and write the date and validity of sterilization on the packaging;

After sterilization, the sterile packaging must be checked for any damage. Sterilization indicators must also be verified. The validity of sterilization must be verified once devices have been sterilized. Do not use the devices if the sterilization period has expired;

During the sterilization process, ensure that the maximum load allowed for the autoclave chamber does not exceed the values recommended in the manufacturer’s instructions;

After sterilization, the devices must only be opened at the moment of use and handled in sterile surgical fields;

It is recommended that the moist heat sterilization method involving the use of an autoclave be used, adhering to the parameters and procedures

established in Plenum's internal protocols and the standard ISO 17665-1. However, the hospital/health service establishment is solely responsible for the methods for cleaning and microbial load control, sterilization equipment, process controls and sterilization procedures used;

Recommended parameters for moist heat sterilization: 121°C for 20 min and 30 min drying time. It is important to note, however, that each user must first read and verify the parameters for their equipment (autoclave) and implement the procedure according to the instruction manual for the respective brand and model of autoclave used by the designated health professional;

Validity of sterilization: 7 to 15 days, provided that the device is packaged in a clean, dry, well-ventilated area away from the action of weather;

ATTENTION: Do not perform autoclaving of this device in its original packaging.

15. EXPIRATION DATE

- Indicated on device labeling.
- The device has an unexpired expiration date.

16. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of Prosthetic Components and comply with health surveillance requirements, the dentist or someone on their team must keep information about the component used in the patient's medical record.

In cases relating to distributors, the information must be passed on to the manufacturer in order to complete the traceability cycle. The information necessary for traceability is related to the device used and surgical procedure as follows:

- Name of the dentist;
- Name of hospital/clinic/office;
- Name of manufacturer and distributor;
- Date of surgery;
- Component code used and batch number;
- Quantity used;
- ANVISA Registration Number.

The responsible dentist and their team should make use of the information contained in the O3 traceability labels supplied with the Prosthetic Component, fixing this information in the patient's clinical record to maintain traceability. In addition, this information must be provided to the patient so that they have information about the device used in their procedure.

The labels contain the following information necessary for the traceability of the prosthetic component used.

- Manufacturer identification;
- Name or commercial model of the device;
- Device code and batch number;
- ANVISA Registration Number.

Traceability information is necessary for notification by the health service and/or by the patient to the National Health Surveillance Agency - ANVISA and to the manufacturer when complaints, technical complaints and/or adverse events occur, for the conduct of the appropriate investments.

17. STORAGE, TRANSPORT, AND HANDLING CONDITIONS

The device packaging must be intact at the moment at which the device is received. Do not use the device if the packaging is damaged or shows signs of tampering. If the packaging is damaged, the device must be discarded;

The data for identification and traceability are on the device label;

The device must be stored and transported in a clean environment, away from sources of heat or humidity, sunlight and the action of the weather in order to prevent any damage or changes in the packaging or its physical form that renders its use impossible;

The device must be stored on shelving or drawers/ cabinet in a manner that allows it to be easily viewed and handled by the operator. Do not store the device near lamps, so that the information contained in the packaging is not affected;

The handling of the device should only be done by qualified and trained personnel to perform this procedure.

18. DISPOSAL OF DEVICE

After use or possible removal of the device from the patient, the prosthetic component that is not subject to further analysis or study must be discarded by the hospital/health service establishment;

It is recommended that the device be completely featureless, so that it cannot be reused. The de-featuring of the device is the sole responsibility of the health service institution /establishment;

After de-featuring, the device must be identified as being unfit for use and disposed of in accordance with health service waste management standards (hospitals and clinics), in accordance with applicable current environmental laws;

The device rendered useless during surgery or laboratory processing must also be discarded following the guidelines described above.

19. CUSTOMER SERVICES

-If the device presents any adverse event with a risk that may harm the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is causing any dissatisfaction, the health professional must generate a technical report containing the details of the non-compliance, with all the claimant data, device code, batch code, and patient data; directly notifying the Plenum Customer Service;

If necessary, send the clean, decontaminated, sterilized, packaged, and properly identified device, with the description of non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - CEP 13212-213 - Jundiaí - SP - Brazil. Phone: +55 11 3109-9000; E-mail: customerservice@plenum.bio

ANVISA also provides a channel to receive notifications of incidents, adverse events, and technical complaints related to the use of devices and services under health surveillance, which can be found on the National Health Surveillance Agency website - ANVISA, at www.anvisa.gov.br, link NOTIVISA.

FDA also provides a channel available to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found in the FDA Manufacturer and User Facility Device Experience (MAUDE) database.

20. LIMITED WARRANTY

Devices manufactured by Plenum are provided with a factory warranty protecting against manufacturing defects. Device defects must be immediately reported to the manufacturer in adherence to legally established deadlines, with an attached technical report detailing the non-compliance identified in the device. This report must be prepared and signed by a qualified health professional and addressed to regulatory authority standards (see the Customer Services section) when it is observed that issues involve a risk to the patient's health.

The warranty of devices manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings, and other information contained in the instructions for use. The lack of such compliance and monitoring will exempt the manufacturer and distributor from any form of liability, particularly in the following cases:

- Reprocessing or reuse of the device by the user;
- Improper use or use other than in a surgical environment.

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