PLENUM[®]
BONE *CRANIUM*[™]
Português

PLENUM[®]
BONE *CRANIUM*[™]
English

PLENUM[®]
BONE *CRANIUM*[™]
Español

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Plenum[®] Bone *cranium*[™]

INSTRUCTIONS FOR USE

TECHNICAL NAME: Bone graft and associated devices

TRADE NAME: Plenum[®] Bone *cranium*[™]

MODEL: For a detailed description of the device, see Table 1.

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E-mail: customerservice@plenum.bio

Website: www.plenum.bio

ATTENTION: Read all the Instructions for Use carefully before use.

Comply with all indications, warnings and precautions mentioned in these Instructions for Use. Failure to observe these points may lead to the occurrence of complications. Product for Medical and Dental Use.

Sterile Product. Sterilized by Gamma Irradiation. Validity: See on the product labeling. Single Use Product. Reprocessing, Resterilizing and Reusing are prohibited.

ANVISA REGISTRATION No.: 81684340019

Code, Batch No., and Date of Manufacture and Sterilization: See the product labeling.

SUMMARY


1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	3
3. DEVICE PRESENTATION FORM	3
4. DEVICE INDICATION, PURPOSE OR INTENDED	4
5. PRINCIPLE OF INTENDED DEVICE FUNCTION AND PERFORMANCE	4
6. INSTRUCTIONS FOR USE	5
7. PRECAUTIONS, RESTRICTIONS AND WARNINGS	6
8. CONTRAINDICATIONS	7
9. POSSIBLE ADVERSE EFFECTS	7
10. SYMBOLS DESCRIPTION	8
11. STERILIZATION	8
12. EXPIRATION DATE	8
13. IDENTIFICATION AND TRACEABILITY	8
14. STORAGE, TRANSPORT, CONSERVATION AND HANDLING CONDITIONS	9
15. DEVICE DISPOSAL	9
16. CUSTOMER SERVICE	9
17. LIMITED WARRANTY	9

1. DETAILED DEVICE DESCRIPTION

The Plenum® Bone *cranium*™ bone graft is a synthetic biphasic bioceramic composed of hydroxyapatite and tricalcium β-phosphate (HA:β-TCP) in a 70:30% ratio. The macro and microstructure, surface area (~2 m²/g) and porosity (~50%) of Plenum® Bone *cranium*™ confer effective osteoconductivity, high hydrophilic property, controlled absorption of the constituent phases of bioceramics, which are absorbed by the body, and replaced by newly formed tissue during the process of bone tissue repair or regeneration. The porosity of Plenum® Bone *cranium*™ ensures efficient osteogenic cell growth by promoting bone regeneration at the site of the bone defect. In turn, the microporosity and surface roughness of Plenum® Bone *cranium*™ facilitate the diffusion of biological fluids and cell adhesion. In addition, the macro and microstructure, surface area, percentage of porosity and the chemical composition of Plenum® Bone *cranium*™ give the product a high hydrophilic property, which makes it easily manipulable (with fluids such as saline solution or the patient’s own blood) in various applications for bone reconstruction. In general, absorbable alloplastic materials for bone repair are the materials of choice; since these materials are free of possible contaminants or organic molecules from homogeneous and heterogenous grafts, which can induce immune responses. In this way, alloplastic materials minimize the risk of infection, inflammation, and other postoperative complications. The Plenum® Bone *cranium*™ bone graft functions as a synthetic matrix for both bulking/rebuilding bone tissue and maintaining the defect space, and is gradually replaced by newly formed bone. The composition of Plenum® Bone *cranium*™ reflects in the controlled absorption behavior. As a faster absorption of β-TCP promotes a continuous increase in the porosity of the material, which allows an integration of the material/tissue due to cell adhesion and proliferation and growth of blood vessels (subsequent deposition of the newly formed bone), HA provides stability of the newly formed bone volume for a longer period.

Plenum® Bone *cranium*™ is a single-use product, supplied in sterile condition and is available in the codes, dimensions, and quantities described in Table 1.

TABLE 1. Illustrative image, codes, dimensions and quantities of the Plenum® Bone *cranium*™

Illustrative image		
Code	Model / Presentation (granules size/quantity)	Model content
018-25-B11-25-05	0.25 - 0.5 mm - 0.5 g	1 unit of 0.5 g
018-25-B11-25-10	0.25 - 0.5 mm - 1.0 g	1 unit of 1.0 g
018-25-B11-25-20	0.25 - 0.5 mm - 2.0 g	1 unit of 2.0 g
018-25-B11-25-40	0.25 - 0.5 mm - 4.0 g	1 unit of 4.0 g
018-25-B11-25-80	0.25 - 0.5 mm - 8.0 g	1 unit of 8.0 g
018-25-B11-50-05	0.5 - 1.0 mm - 0.5 g	1 unit of 0.5 g
018-25-B11-50-10	0.5 - 1.0 mm - 1.0 g	1 unit of 1.0 g
018-25-B11-50-20	0.5 - 1.0 mm - 2.0 g	1 unit of 2.0 g
018-25-B11-50-40	0.5 - 1.0 mm - 4.0 g	1 unit of 4.0 g

018-25-B11-50-80	0.5 - 1.0 mm - 8.0 g	1 unit of 8.0 g
018-25-B11-10-05	1.0 - 2.0 mm - 0.5 g	1 unit of 0.5 g
018-25-B11-10-10	1.0 - 2.0 mm - 1.0 g	1 unit of 1.0 g
018-25-B11-10-20	1.0 - 2.0 mm - 2.0 g	1 unit of 2.0 g
018-25-B11-10-40	1.0 - 2.0 mm - 4.0 g	1 unit of 4.0 g
018-25-B11-10-80	1.0 - 2.0 mm - 8.0 g	1 unit of 8.0 g
1350705	NewCovery Bone 0.25 - 0.5 mm - 0.5 g	1 unit of 0.5 g
1350706	NewCovery Bone 0.25 - 0.5 mm - 1 g	1 unit of 1.0 g
1350707	NewCovery Bone 0.25 - 0.5 mm - 2 g	1 unit of 2.0 g
1350715	NewCovery Bone 0.25 - 0.5 mm - 4 g	1 unit of 4.0 g
1350718	NewCovery Bone 0.25 - 0.5 mm - 8 g	1 unit of 8.0 g
1350708	NewCovery Bone 0.5 - 1 mm - 0.5 g	1 unit of 0.5 g
1350709	NewCovery Bone 0.5 - 1 mm - 1 g	1 unit of 1.0 g
1350710	NewCovery Bone 0.5 - 1 mm - 2 g	1 unit of 2.0 g
1350716	NewCovery Bone 0.5 - 1 mm - 4 g	1 unit of 4.0 g
1350720	NewCovery Bone 0.5 - 1 mm - 8 g	1 unit of 8.0 g
1350711	NewCovery Bone 1 - 2 mm - 0.5 g	1 unit of 0.5 g
1350712	NewCovery Bone 1 - 2 mm - 1g	1 unit of 1.0 g
1350714	NewCovery Bone 1 - 2 mm - 2 g	1 unit of 2.0 g
1350717	NewCovery Bone 1 - 2 mm - 4 g	1 unit of 4.0 g
1350721	NewCovery Bone 1 - 2 mm - 8 g	1 unit of 8.0 g

2. DEVICE COMPOSITION

Plenum® Bone *cranium*™ is an absorbable synthetic bone substitute based on biphasic bioceramics composed of hydroxyapatite and tricalcium β-phosphate (HA:β-TCP) in a 70:30% ratio.

3. DEVICE PRESENTATION FORM

The Plenum® Bone *cranium*™ has no accessories.

Plenum® Bone *cranium*™ is supplied packaged in the models described in the table above, in sterile condition, and presented as follows:

0.5, 1.0 AND 2.0 G MODELS:

Primary packaging: polyethylene terephthalate (PET) blister and thermo-sealed with Tyvek-type surgical grade paper;

Secondary packaging: case with transparent lid in crystal polystyrene (GPPS) and white base in polypropylene (PP), containing 01 blister of the product in one of the presentations described above;

Tertiary packaging: in the form of a PET blister and thermo-sealed with surgical grade Tyvek paper, sterilized by Gamma Irradiation;

Outer packaging: rigid carton paper envelope containing primary, secondary, tertiary packaging and 3 identification and traceability labels;

4.0 AND 8.0 G MODELS:

Primary packaging: polyethylene terephthalate (PET) blister and heat-sealed with surgical grade Tyvek paper;

Secondary packaging: PET blister and thermo-sealed with Tyvek-type surgical grade paper, sterilized by Gamma Irradiation;

Outer packaging: rigid carton paper envelope containing primary, secondary, tertiary packaging and 3 identification and traceability labels;

It is the responsibility of the hospital institution/health service establishment to ensure the information for the traceability of the Plenum® Bone *cranium*™, by annotating in the patient's medical record the code and batch of the product used, or by means of one of the 3 traceability labels that accompany the package;

If you want to obtain the Instructions for Use in printed format and at no additional cost (including shipping), request it by calling (11) 3109-9000, or through the electronic address (relacionamento@plenum.bio) or through the website (www.plenum.bio/ifu);

See below the illustrative images of the primary, secondary and/or tertiary (Figure 1) and external (Figure 2) packaging of a Plenum® Bone *cranium*™ model in the form in which it is marketed.

FOR 0.5, 1.0 AND 2.0 G MODELS:

FOR 4.0 AND 8.0 G MODELS:


FIGURE 1: ILLUSTRATIVE IMAGES OF PRIMARY, SECONDARY AND/OR TERTIARY PACKAGING.

FIGURE 2: ILLUSTRATIVE IMAGES OF THE OUTER PACKAGING.

4. DEVICE INDICATION, PURPOSE OR INTENDED

The Plenum® Bone *cranium*™ is a synthetic bone graft, which is intended for the filling and/or reconstruction of bone defects in the cranio-maxillofacial region, and is indicated for guided tissue regeneration, alveolar ridge augmentation/reconstruction, maxillary sinus floor lifting or bone defects resulting from congenital, post-traumatic, post-surgical problems, or involvement of bone structures due to diseases such as tumors or infectious processes, which are not intrinsic to the stability of the bone structure.

This synthetic bone graft works as a bone substitute, favoring the maintenance of the bone defect space, providing the repair of the defect or even the increase in bone volume required by the surgeon. The Plenum® Bone *cranium*™ is absorbed by the body, and replaced by newly formed tissue during the process of bone tissue repair or regeneration.

5. PRINCIPLE OF INTENDED DEVICE FUNCTION AND PERFORMANCE

The Plenum® Bone *cranium*™ has a macro, micro and nanostructure that mimics the trabecular bone tissue providing surface area, percentage of porosity and unique physicochemical characteristics that promote excellent osteoconductive property, high hydrophilicity and adaptability to the recipient bed. Characteristics that facilitate the diffusion of biological fluids and cell adhesion to the product, being ideal for large reconstructions that require the bone graft to promote vascularization and nutrition in all its scope. The absorption of Plenum® Bone *cranium*™ occurs in a controlled manner of the constituent phases of bioceramics, which are absorbed by the body and replaced by newly formed tissue during the process of repair or regeneration of bone tissue.

The intended performance of Plenum® Bone *cranium*™ is that it functions as an absorbable synthetic bone substitute for the filling and/or reconstruction of bone defects in the craniomaxillofacial area, providing space maintenance and repair/regeneration of bone defects.

6. INSTRUCTIONS FOR USE

Carefully read all instructions for use, recommendations and guidelines before using the Plenum® Bone *cranium*™ synthetic bone graft;

Before use, check the expiration date of the product. Do not use with the packaging violated or with the expiration date;

The Plenum® Bone *cranium*™ is supplied ready for use. It is available in different quantities and in three particle size ranges (0.25 – 0.50 mm, 0.50 – 1 mm and 1 -2 mm). The surgeon can perform prior hydration with saline solution or the patient’s own blood for effectiveness in handling the product;

Plenum® Bone *cranium*™ is available in various models and sizes of granules, and it is the surgeon’s responsibility (doctor or dentist) to choose the product model with the most appropriate characteristics and dimensions, after the individual clinical analysis of each case;

The techniques for using and applying Plenum® Bone *cranium*™ vary according to the surgeon’s preference (physician or dentist), and it is up to him to choose the therapeutic approach, model and dimensions of the product, application technique, use of complementary biomaterials, as well as the criteria for monitoring and evaluating the results of the surgery;

Proper preparation of the site (recipient bed) where the Plenum® Bone *cranium*™ will be applied is very important. The entire area to be treated must be filled with the product, as well as bleeding must be controlled, so that the bone graft has a perfect adaptation to the recipient bed and achieves the expected results;

In the procedure of manipulation and implantation of the Plenum® Bone *cranium*™, the surgeon should use only appropriate and sterile surgical instruments;

Before the surgical procedure, the surgeon (doctor or dentist) must carry out a rigorous preoperative planning, through clinical examinations and images (radiographs or CT scans). You should also consider the need for a detailed anamnesis and complementary tests on the patient’s general health (complete blood count, coagulogram, calcium dosage, etc.), since this information directly influences the body’s biological response to the use of the bone graft, Plenum® Bone *cranium*™, and then previously define a surgical plan, with the selection of the size of the granules and the most appropriate quantity for the patient’s condition;

The surgical technique to be used should be indicated by the surgeon, who should consider his clinical and surgical experience and choose the most appropriate technique for each patient (individualized procedure). Therefore, a general surgical technique is not recommended for all patients.

The surgical protocol should be carried out according to the surgeon’s previous references and experiences, always considering the most appropriate choice of granule size and quantity of the product, technique, grafting sequence, and the use of membranes, meshes, or meshes when necessary, always based on conventional and established therapeutic techniques of tissue regeneration and bone grafting;

Postoperative care and correct hygiene/cleanliness of the operated site are extremely important. The patient should be instructed about the existing risks until there is complete repair of the bone and/or gingival tissue at the surgical site;

The patient should be adequately oriented regarding postoperative care. The patient’s ability and willingness to follow the surgeon’s instructions is one of the most important aspects in a surgical procedure;

For elderly patients, those with mental problems, drug addicts, or with any other difficulty in following postoperative recommendations, this product may represent a higher risk of failure, as these patients may ignore the instructions and restrictions arising from this procedure;

The surgeon should adopt a correct clinical/surgical planning described in the specific literature on bone regeneration and bone grafting techniques for the immediate and late pre- and postoperative periods;

The surgeon must clearly guide the patient, who must understand the importance of follow-up in the postoperative period, in which the stability and integrity of the product must be verified, as well as the incorporation and/or absorption at the appropriate time;

In order to protect the patient and his surgeon, it is important to ensure the correct and total traceability of the Plenum® Bone *cranium*™ used, by annotating the product code and lot in the patient protocol, or by fixing the identification and traceability labels that accompany the product in the patient protocol.

HANDLING AND MANIPULATION OF THE PLENUM® BONE CRANIUM™


0.5, 1.0 AND 2.0G MODELS:

Plenum® Bone *cranium*™ product in the form of presentation in the primary and secondary packaging.


Turn the lid counterclockwise as indicated in the OPEN direction on the lid.


Remove the lid to remove the product from its primary packaging. Place the cap over the sterile drape for use with the Plenum® Bone *cranium*™ product handling container.

Open the primary pack (blister) with your thumb and index fingers holding the flap of the surgical paper. Be careful with sudden movements so as not to drop the product at the time of opening.


Pour the Plenum® Bone *cranium*™ product into the handling container, which can be the lid of the sterile case (secondary packaging), to perform previous hydration steps with saline solution or the patient’s own blood.

Examples of bone defect reconstruction in the craniomaxillofacial region, maxillary sinus and inferior alveolar ridge, filled with the Plenum® Bone *cranium*™ product for all models.

ATTENTION: The figures are merely illustrative. They do not represent the real dimensions and characteristics of the product.


4.0 E 8.0 G MODELS:

Plenum® Bone *cranium*™ product in the form of presentation in the primary packaging.

Open the primary pack (blister) with your thumb and index fingers holding the flap of the surgical paper. Be careful with sudden movements so as not to drop the product at the time of opening.

Perform previous hydration steps with saline solution or the patient's own blood in the primary packaging of the Plenum® Bone *cranium*™.

The Plenum® Bone *cranium*™ can be manipulated in its own packaging.

ATTENTION: The figures are merely illustrative. They do not represent the real dimensions and characteristics of the product.

7. PRECAUTIONS, RESTRICTIONS AND WARNINGS

The surgeon should not use Plenum® Bone *cranium*™ until he has carefully read all instructions for use and the information contained in the package. All usage information, warnings, precautions, and recommendations referenced in these instructions for use should be considered. Failure to observe these points may lead to the occurrence of complications for the surgeon and especially for the patient;

The handling of the Plenum® Bone *cranium*™ should be performed exclusively by surgeons (physicians or dentists) previously trained and qualified for this procedure, with experience and knowledge about guided tissue regeneration techniques and the use of bone substitutes;

It is recommended to request preoperative tests according to the patient's anamnesis;

Plenum® Bone *cranium*™ is supplied in a sterile condition, sterilized by Gamma Irradiation. Do not use the product if the packaging is violated or expired;

SINGLE-USE product, therefore, reprocessing, resterilization and reuse are prohibited and may bring risks of failure and infection to the patient. If the product is exposed to one of these situations, it will lead to contamination and permanent damage to the product and/or cross-infection in the patient;

The product label presents the expiration date, the code and batch number of the product, ANVISA registration number, technical responsible and the company's data, as well as the respective symbology, allowing full traceability of the product;

The product should be kept in the original packaging until the moment of use, and should only be handled in a sterile surgical drape;

Any product that is contaminated or damaged during handling and preparation must be discarded according to the guidelines described in the item Product Disposal;

The surgeon must understand and plan all the steps of the procedure and verify the amount he intends to use. Always have the product available in quantities higher than the one you intend to use;

The surgeon must be aware of the patient's clinical conditions. In the preoperative examination, the site where the Plenum® Bone *cranium*™ will be placed must be rigorously verified through clinical and imaging examinations. The surgeon must consider the general clinical status of the patient, analyze local or systemic factors that may interfere with the process of bone healing, soft tissues, and incorporation and absorption of bone substitute;

Before starting surgery, the surgeon must submit the patient to a rigorous oral antisepsis to prevent the product from coming into contact with contaminants;

The patient should be warned that failure to comply with the instructions of the postoperative period may lead to the failure of the surgical procedure, requiring new surgery for revision or removal of the bone substitute;

Patients with inadequate oral hygiene, clinical conditions with hematological disorders and other uncontrolled systemic alterations, have restrictions on the use of this product;

Forbidden to use any leftovers of Plenum® Bone *cranium*™. The product that may be left over after the end of a surgery must be disposed of in accordance with the waste management standards of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation in force;

The wrong choice of the bone graft that will be used, as well as errors in the indication, manipulation, preparation of bone tissue or soft tissues, and

installation of the product, can cause damage to the physical structure and contamination of the product and contribute to its failure;

The use of this bone graft with inappropriate techniques, procedures, and surgical conditions may harm the patient, causing unsatisfactory results in tissue repair or regeneration;

The main and most probable causes of the lack of incorporation and absorption of the Plenum® Bone *cranium*™ with the consequent loss of surgery are: infections, localized or systemic diseases, inappropriate use of membranes, meshes or meshes, insufficient blood or nutrient supply, inadequate hygiene, trauma or overload at the surgical site, contamination of the product in the pre and/or intraoperative or postoperative period when there is exposure of the bone substitute in the oral cavity, mechanical instability of the granules, bone necrosis, lack of professional training, technical error in the installation of the product, and patients who do not follow the postoperative instructions. The surgeon should consider all these aspects in his preoperative planning;

Plenum® Bone *cranium*™ has not been tested in lactating women, therefore it is not known whether the metabolites of Plenum® Bone *cranium*™ pass into breast milk;

Plenum® Bone *cranium*™ is only an auxiliary material for bone replacement. Use alone does not guarantee the success of the treatment. The success of the therapy depends on the intervention technique and the patient's reaction;

The success of a bone reconstruction/regeneration surgical procedure using Plenum® Bone *cranium*™ associated with other substitutes, bone grafts or membranes is related to the correct indication and selection of the product, the correct surgical technique, the preparation, the stability of the product, the time and quality of the biomaterial(s)/recipient bed incorporation, as these extrinsic factors can be affected by numerous biological factors, clinical and biomechanical. These factors can cause failure and early loss of the product and compromise the surgical procedure. Therefore, faithful obedience to all indications for use, contraindications, precautions, warnings and other observations contained in these instructions for use are essential for the success of the surgical procedure;

The patient should be informed of all the limitations, risks, complications, warnings, and possible adverse effects arising from the guided tissue regeneration procedure and bone grafting techniques with the use of synthetic bone substitutes;

The patient should be informed of the risks involved in procedures with the use of bone substitutes where the operated site will suffer stress and direct mechanical stress, for example, exposure to constant deformation forces and/or tension;

The patient should be informed of the risks involved in the application of excessive mechanical load or trauma affecting the operated region, which may imply failure, loss of mechanical stability, exposure, and displacement of the particles of the synthetic bone substitute;

It is important and fundamental that the surgeon explains all the above topics to his patient, and that he has a complete and clear understanding of all these aspects.

8. CONTRAINDICATIONS

The surgeon should evaluate the patient's condition prudently before surgery and the use of Plenum® Bone *cranium*™, and check for any contraindications listed below. Contraindications may include, but are not limited to:

Patient who has insufficient bone, gingival or soft tissue quantity and quality to provide primary stability of the synthetic bone graft;

Patient with poor bone quality due to osteoporosis or osteopenia, which

can cause early failure of the bone graft due to inefficient tissue repair or regeneration;

Use in places that suffer stress, direct mechanical stress, or the granules are exposed to constant deformation forces and/or tension;

Patients who do not have sufficient soft tissues for complete coverage of the bone graft;

Inflammation and/or infection near or at the surgical site;

Recent history of systemic or localized infection;

Fever;

Leukocytosis;

Patient with allergy or sensitivity to the raw material used in the manufacture of Plenum® Bone *cranium*™ (see "PRODUCT COMPOSITION");

Patient with developmental disorders of the skeletal system;

Pregnancy and lactation period;

Use in patients with immunological alterations and systemic diseases that result in poor healing or the absence of integrity of the existing tissue at the surgical site;

Patient with vascular disorders or coagulation disorders;

Calcium metabolism disorder;

Treatment with steroids and other drugs that interfere with calcium metabolism;

Alcoholic patient or drug user who cannot correctly follow postoperative instructions;

Patient with mental illness or unable to follow the recommendations and care of the postoperative phase;

Patient unable to maintain good hygiene at the surgical site;

Immunosuppressive pathologies;

Patient submitted to radiotherapy or use of drugs that compromise the vascularization of tissues, bone and soft;

All known contraindications in the medical-scientific literature for the use of absorbable synthetic bone substitutes.

9. POSSIBLE ADVERSE EFFECTS

Positive results may not be achieved in all surgical procedures. Every surgical procedure has risks and possibilities of complications, such as infections, hemorrhages, allergic drug reactions and anesthetic risks, among others. Possible adverse effects associated with the use of Plenum® Bone *cranium*™ may include, but are not limited to:

The surgical procedure for the installation of bone grafts or substitutes can bring complications in the trans and postoperative periods, such as: hemorrhage, pain, edema, abscess, fistula, dehiscence, nerve injury or compression (sensory or motor paresthesia, partial or total, transient or permanent), soft tissue injury, local or systemic infection, and other complications typical of a procedure for the placement of bone grafts or substitutes;

Sensitivity of the tissue at the surgical site;

Pain or discomfort due to the presence of the bone graft/substitute or due to the surgical procedure;

Lack of incorporation of the bone graft, or even delay in healing, which can result in product failure;

Loss of stability, displacement of the granules of the bone graft resulting from non-incorporation and early partial absorption, or even caused by inadequate installation of the product;

Allergic reaction to the raw material of Plenum® Bone *cranium*™;

Intense fibrotic tissue reaction around the surgical area;

Bone necrosis, insufficient vascularization, and insufficient bone formation, which can cause loss of the grafted area;

Premature dehiscence of sutures, tension in the coaptation of the surgical edges, or premature exposure of the Plenum® Bone *cranium*™ and the grafted area can alter scarring processes and result in failure of the surgical procedure;


The patient must be informed and understand that in all surgical procedures there is the possibility of complications. The patient should be informed by the surgeon of all these complications, and of the risks arising from tissue regeneration surgeries and bone grafting techniques with the use of absorbable synthetic bone grafts/substitutes.

10. SYMBOLS DESCRIPTION

Table 2 describes the symbols that may be on product labeling. Please check the physical labeling for the symbols applied to the product. The symbols express information related to the identification and use of the product.

TABLE 2: Symbols used in labeling to express information related to the identification and use of the product.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using Irradiation
	Do not resterilize
	Do not use if package is damaged
	Double sterile barrier system
	Keep away from sunlight
	Keep dry
	Do not re-use
	Consult electronic instructions for use (plenum.bio/ifu)
	Unique device identifier
	Device manufacturing material

	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
	Reprocessing prohibited (Brazilian Requirement)
	Consult the instructions for use in electronic format (QR code)
	ANVISA registration number

NOTE: The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

11. STERILIZATION

Plenum® Bone *cranium*™ is supplied in sterile condition, sterilized by Gamma Irradiation;

Single-use medical device and cannot be reprocessed, re-sterilized or reused;

Before use, check the integrity of the package and the expiration date. Do not use the product if the packaging is violated or the product has expired;

Date of manufacture, expiration date, product code, batch number and other information regarding identification and traceability are on the product labeling;

To maintain sterility, the product should be opened only at the time of use and handled in appropriate and aseptic environments.

The opening of the package for the use of the product should only be done by qualified and trained personnel to perform this procedure.

12. EXPIRATION DATE

Check the expiration date of the product on the label.

The date of manufacture, expiration date, code, and batch number of the product must be consulted on the product labeling;

Before use check the validity, do not use the product if it has expired.

13. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of the Plenum® Bone *cranium*™ and comply with the requirements of current legislation, the surgeon or someone from his team must keep in the patient’s medical record the information regarding the product used. In addition, this information must also be passed on to the product distributor and the patient, in order to complete the traceability cycle. The information required for traceability is that relating to the product used and surgery, as follows:

- Name of the surgeon;
- Name of the hospital/clinic/office;
- Name of the manufacturer and distributor;
- Date of surgery;
- Code of the product used;
- Product batch number;

- Quantity used;
- Product registration with ANVISA.

The surgeon in charge and his team must make use of the information contained in the 3 traceability labels provided with the Plenum® Bone *cranium*™, fixing this information in the patient's clinical record for the maintenance of traceability. In addition, this information must be provided to the patient so that he has information about the product used in his procedure.

The labels contain the following information necessary for the traceability of the product used:

- Manufacturer identification;
- Name or business model of the product;
- Product code;
- Product batch number;
- ANVISA registration number.

Traceability information is necessary for notification by the health service and/or by the patient himself to the National Health Surveillance Agency – ANVISA and to the manufacturer, when adverse events occur, in order to conduct the appropriate investigations.

14. STORAGE, TRANSPORT, CONSERVATION AND HANDLING CONDITIONS

The packaging must be intact at the time of receipt, do not use the product if the packaging is damaged or violated, or with the expiration date expired. If the packaging is damaged, the product should be considered as non-sterile and discarded;

The product must be stored and transported in a clean environment (maximum temperature limit of up to 45 °C), protected from sources of heat or humidity, away from the incidence of sunlight and free from the action of bad weather, in order to prevent any damage or alteration in the packaging or its physical form that makes it impossible to use;

The storage of the product on shelves or drawers/cabinet should be in such a way that it facilitates the operator's viewing and handling. Do not store the product near lamps, so that the information contained in the package is not affected;

The product should be opened only at the time of use in a sterile surgical field;

The handling of the product should only be done by qualified and trained personnel to perform this procedure.

15. DEVICE DISPOSAL

Discard the product if any of the successive packages that ensure sterility are damaged or opened;

After the use or eventual removal of the product from the patient, the synthetic bone graft (Plenum® Bone *cranium*™) that is not the object of further analysis or studies must be discarded by health service establishments (hospitals, offices and clinics);

It is recommended that the Plenum® Bone *cranium*™ be completely de-characterized, so that it cannot be reused. The decharacterization of the product is the sole responsibility of the hospital institution/health service establishment;

After decharacterization, the Plenum® Bone *cranium*™ must be identified as being unfit for use and disposed of in accordance with the standards of waste management of health services (hospitals, offices and clinics) in accordance with applicable environmental legislation;

Plenum® Bone *cranium*™ that is left over and not used during surgery should also be discarded following the guidelines described above.

16. CUSTOMER SERVICE

If the product presents any adverse event, with potential risk that may injure the patient, generates or has the potential for injury or threat to public health, is out of its specifications, or is generating any dissatisfaction, the health professional must generate a technical report containing the details of the non-conformity, with all the claimant's data, product code and lot, and patient data; by directly notifying Plenum's Customer Relationship Service;

If necessary, send the clean, decontaminated and packaged product, duly identified, and with a description of the non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lote 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logístico - CEP 13212-213 - Jundiaí - SP - Brazil. Phone: (11) 3109-9000; Email: relacionamento@plenum.bio

ANVISA also provides a channel to receive notifications of incidents, adverse events and technical complaints, related to the use of products and services under health surveillance, which can be found on the website of the National Health Surveillance Agency – ANVISA, at the electronic address www.anvisa.gov.br, link NOTIVISA.

17. LIMITED WARRANTY

The products manufactured by Plenum have a factory warranty against manufacturing defects. The occurrence of defects in the product must be immediately reported to the manufacturer, observing the legal deadline, attaching a technical report with the details of the non-conformity found. This report must be generated and signed by a qualified and trained health professional, also observing the relevant ANVISA rules in case of risk to the patient's health (see the Customer Service topic above).

The Warranty of products manufactured by Plenum is subject to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings and other information contained in the instructions for use. Failure to comply with this and follow exempts the manufacturer and distributor from any type of liability, especially in the following cases:

- Reprocessing, re-sterilization or reuse of the product by the user;
- Inappropriate use or other use other than in a surgical environment;
- If the product has not been used in full compliance with the instructions for use.

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