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TEMPORARY
ABUTMENT
English

PILAR
PROVISIONAL
Español

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TEMPORARY ABUTMENT

INSTRUCTIONS FOR USE

TECHNICAL NAME: Dental Implant Components

TRADE NAME: Temporary Abutment

MODEL: For a detailed description of the device, see Table 1.

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Website: www.plenum.bio

ATTENTION:

Carefully read these instructions for use before using the device. Follow the recommendations, warnings, and precautions mentioned in these Instructions for Use. A failure to observe these items may result in complications.

Device for Medical and Dental use. Sterile device. Sterilized using Gamma Irradiation. Expiration date: See device labeling. Single-use only. Prohibited reprocessing, resterilize, and reuse.

ANVISA DEVICE REGISTRATION N°.: 81684340012

Device code, batch code, and date of manufacture: See device labeling.

REV. 08

SUMMARY

1. DETAILED DEVICE DESCRIPTION	3
2. DEVICE COMPOSITION	6
3. DEVICE PRESENTATION FORM	6
4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE	6
5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE	7
6. INSTRUCTIONS FOR USE	7
7. PRECAUTIONS	8
8. WARNINGS	9
9. CONTRAINDICATIONS	9
10. MAGNETIC RESONANCE IMAGING	9
11. POSSIBLE ADVERSE EFFECTS	10
12. SYMBOLS DESCRIPTION	10
13. AUXILIARY INSTRUMENTS	10
14. STERILIZATION	11
15. EXPIRATION DATE	11
16. IDENTIFICATION AND TRACEABILITY	11
17. STORAGE, TRANSPORT AND HANDLING CONDITIONS	11
18. DISPOSAL OF DEVICE	11
19. CUSTOMER SERVICES	12
20. LIMITED WARRANTY	12

1. DETAILED DEVICE DESCRIPTION

The Plenum’s Temporary Abutment line of devices is intended to be installed as an intermediate prosthetic component between the implant and the temporary dental prosthesis while awaiting tissue repair after the dental implant has been installed or a permanent prosthesis is being prepared.

The devices made available in the models described in the table below must only be used with Plenum’s Dental Implants and Instruments (Implants and Instruments are not described in these instructions for use and must be purchased separately).

Illustrative images, codes, models, prosthetic platforms, and main dimensions of Temporary Abutments are presented in Table 1.

TABLE 1: Illustrative images, codes, models, prosthetic platforms, and dimensions of Temporary Abutments.

Illustrative Image	Codes*	Models	Model Composition	Diameter (mm)	Length (mm)	Prosthetic Interface
	018-21-C60-33-00	SL Ø3.5X0MM	SL Ø3.5 X 0MM + HEXAGONAL SCREW SL 0MM	3.5	0.0	SL
	018-21-C60-33-01	SL Ø3.5X1MM	SL Ø3.5 X 1MM + HEXAGONAL SCREW SL 1MM	3.5	1.0	SL
	018-21-C60-33-02	SL Ø3.5X2MM	SL Ø3.5 X 2MM + HEXAGONAL SCREW SL 2MM	3.5	2.0	SL
	018-21-C60-33-03	SL Ø3.5X3MM	SL Ø3.5 X 3MM + HEXAGONAL SCREW SL 3MM	3.5	3.0	SL
	018-21-C60-33-04	SL Ø3.5X4MM	SL Ø3.5 X 4MM + HEXAGONAL SCREW SL 4MM	3.5	4.0	SL
	018-21-C60-33-05	SL Ø3.5X5MM	SL Ø3.5 X 5MM + HEXAGONAL SCREW SL 5MM	3.5	5.0	SL
	018-21-C60-33-06	SL Ø3.5X6MM	SL Ø3.5 X 6MM + HEXAGONAL SCREW SL 6MM	3.5	6.0	SL
	018-22-C60-33-00	RE Ø3.6X0MM	RE Ø3.6 X 0MM + HEXAGONAL SCREW RE 0MM	3.6	0.0	RE
	018-22-C60-33-01	RE Ø3.6X1MM	RE Ø3.6 X 1MM + HEXAGONAL SCREW RE 1MM	3.6	1.0	RE
	018-22-C60-33-02	RE Ø3.6X2MM	RE Ø3.6 X 2MM + HEXAGONAL SCREW RE 2MM	3.6	2.0	RE
	018-22-C60-33-03	RE Ø3.6X3MM	RE Ø3.6 X 3MM + HEXAGONAL SCREW RE 3MM	3.6	3.0	RE
	018-22-C60-33-04	RE Ø3.6X4MM	RE Ø3.6 X 4MM + HEXAGONAL SCREW RE 4MM	3.6	4.0	RE
	018-22-C60-33-05	RE Ø3.6X5MM	RE Ø3.6 X 5MM + HEXAGONAL SCREW RE 5MM	3.6	5.0	RE
	018-22-C60-33-06	RE Ø3.6X6MM	RE Ø3.6 X 6MM + HEXAGONAL SCREW RE 6MM	3.6	6.0	RE
	018-22-C60-45-00	RE Ø4.5X0MM	RE Ø4.5 X 0MM + HEXAGONAL SCREW RE 0MM	3.6	0.0	RE
	018-22-C60-45-01	RE Ø4.5X1MM	RE Ø4.5 X 1MM + HEXAGONAL SCREW RE 1MM	3.6	1.0	RE
	018-22-C60-45-02	RE Ø4.5X2MM	RE Ø4.5 X 2MM + HEXAGONAL SCREW RE 2MM	3.6	2.0	RE
	018-22-C60-45-03	RE Ø4.5X3MM	RE Ø4.5 X 3MM + HEXAGONAL SCREW RE 3MM	3.6	3.0	RE
	018-22-C60-45-04	RE Ø4.5X4MM	RE Ø4.5 X 4MM + HEXAGONAL SCREW RE 4MM	3.6	4.0	RE
	018-22-C60-45-05	RE Ø4.5X5MM	RE Ø4.5 X 5MM + HEXAGONAL SCREW RE 5MM	3.6	5.0	RE
	018-22-C60-45-06	RE Ø4.5X6MM	RE Ø4.5 X 6MM + HEXAGONAL SCREW RE 6MM	3.6	6.0	RE
	018-24-C60-33-01	UF Ø3.6X1MM	UF Ø3.6 X 1MM + HEXAGONAL SCREW UF 1MM	3.6	1.0	UF
	018-24-C60-33-02	UF Ø3.6X2MM	UF Ø3.6 X 2MM + HEXAGONAL SCREW UF 2MM	3.6	2.0	UF
	018-24-C60-33-03	UF Ø3.6X3MM	UF Ø3.6 X 3MM + HEXAGONAL SCREW UF 3MM	3.6	3.0	UF
	018-24-C60-33-04	UF Ø3.6X4MM	UF Ø3.6 X 4MM + HEXAGONAL SCREW UF 4MM	3.6	4.0	UF
	018-24-C60-33-05	UF Ø3.6X5MM	UF Ø3.6 X 5MM + HEXAGONAL SCREW UF 5MM	3.6	5.0	UF
	018-24-C60-33-06	UF Ø3.6X6MM	UF Ø3.6 X 6MM + HEXAGONAL SCREW UF 6MM	3.6	6.0	UF
	018-24-C60-45-01	UF Ø4.5X1MM	UF Ø4.5 X 1MM + HEXAGONAL SCREW UF 1MM	4.5	1.0	UF
	018-24-C60-45-02	UF Ø4.5X2MM	UF Ø4.5 X 2MM + HEXAGONAL SCREW UF 2MM	4.5	2.0	UF
	018-24-C60-45-03	UF Ø4.5X3MM	UF Ø4.5 X 3MM + HEXAGONAL SCREW UF 3MM	4.5	3.0	UF
	018-24-C60-45-04	UF Ø4.5X4MM	UF Ø4.5 X 4MM + HEXAGONAL SCREW UF 4MM	4.5	4.0	UF
	018-24-C60-45-05	UF Ø4.5X5MM	UF Ø4.5 X 5MM + HEXAGONAL SCREW UF 5MM	4.5	5.0	UF
	018-24-C60-45-06	UF Ø4.5X6MM	UF Ø4.5 X 6MM + HEXAGONAL SCREW UF 6MM	4.5	6.0	UF

	018-25-C60-33-01	MF Ø3.6X1MM	MF Ø3.6 X 1MM + HEXAGONAL SCREW MF 1MM	3.6	1.0	MF
	018-25-C60-33-02	MF Ø3.6X2MM	MF Ø3.6 X 2MM + HEXAGONAL SCREW MF 2MM	3.6	2.0	MF
	018-25-C60-33-03	MF Ø3.6X3MM	MF Ø3.6 X 3MM + HEXAGONAL SCREW MF 3MM	3.6	3.0	MF
	018-25-C60-33-04	MF Ø3.6X4MM	MF Ø3.6 X 4MM + HEXAGONAL SCREW MF 4MM	3.6	4.0	MF
	018-25-C60-33-05	MF Ø3.6X5MM	MF Ø3.6 X 5MM + HEXAGONAL SCREW MF 5MM	3.6	5.0	MF
	018-25-C60-33-06	MF Ø3.6X6MM	MF Ø3.6 X 6MM + HEXAGONAL SCREW MF 6MM	3.6	6.0	MF
	018-25-C60-45-01	MF Ø4.5X1MM	MF Ø4.5 X 1MM + HEXAGONAL SCREW MF 1MM	4.5	1.0	MF
	018-25-C60-45-02	MF Ø4.5X2MM	MF Ø4.5 X 2MM + HEXAGONAL SCREW MF 2MM	4.5	2.0	MF
	018-25-C60-45-03	MF Ø4.5X3MM	MF Ø4.5 X 3MM + HEXAGONAL SCREW MF 3MM	4.5	3.0	MF
	018-25-C60-45-04	MF Ø4.5X4MM	MF Ø4.5 X 4MM + HEXAGONAL SCREW MF 4MM	4.5	4.0	MF
	018-25-C60-45-05	MF Ø4.5X5MM	MF Ø4.5 X 5MM + HEXAGONAL SCREW MF 5MM	4.5	5.0	MF
	018-25-C60-45-06	MF Ø4.5X6MM	MF Ø4.5 X 6MM + HEXAGONAL SCREW MF 6MM	4.5	6.0	MF
	018-23-C60-33-00	SH Ø4.5X0MM	SH Ø4.5 X 0MM + HEXAGONAL SCREW SH 0MM	4.5	0.0	SH
	018-23-C60-33-01	SH Ø4.5X1MM	SH Ø4.5 X 1MM + HEXAGONAL SCREW SH 1MM	4.5	1.0	SH
	018-23-C60-33-02	SH Ø4.5X2MM	SH Ø4.5 X 2MM + HEXAGONAL SCREW SH 2MM	4.5	2.0	SH
	018-23-C60-33-03	SH Ø4.5X3MM	SH Ø4.5 X 3MM + HEXAGONAL SCREW SH 3MM	4.5	3.0	SH
	018-23-C60-33-04	SH Ø4.5X4MM	SH Ø4.5 X 4MM + HEXAGONAL SCREW SH 4MM	4.5	4.0	SH
	018-23-C60-33-05	SH Ø4.5X5MM	SH Ø4.5 X 5MM + HEXAGONAL SCREW SH 5MM	4.5	5.0	SH
	018-23-C60-33-06	SH Ø4.5X6MM	SH Ø4.5 X 6MM + HEXAGONAL SCREW SH 6MM	4.5	6.0	SH
	018-21-C62-33-00	SMART SL Ø3.5X0MM	SMART SL Ø3.5 X 0MM + HEXALOBULAR SCREW SL 0MM	3.5	0.0	SL
	018-21-C62-33-01	SMART SL Ø3.5X1MM	SMART SL Ø3.5 X 1MM + HEXALOBULAR SCREW SL 1MM	3.5	1.0	SL
	018-21-C62-33-02	SMART SL Ø3.5X2MM	SMART SL Ø3.5 X 2MM + HEXALOBULAR SCREW SL 2MM	3.5	2.0	SL
	018-21-C62-33-03	SMART SL Ø3.5X3MM	SMART SL Ø3.5 X 3MM + HEXALOBULAR SCREW SL 3MM	3.5	3.0	SL
	018-21-C62-33-04	SMART SL Ø3.5X4MM	SMART SL Ø3.5 X 4MM + HEXALOBULAR SCREW SL 4MM	3.5	4.0	SL
	018-21-C62-33-05	SMART SL Ø3.5X5MM	SMART SL Ø3.5 X 5MM + HEXALOBULAR SCREW SL 5MM	3.5	5.0	SL
	018-21-C62-33-06	SMART SL Ø3.5X6MM	SMART SL Ø3.5 X 6MM + HEXALOBULAR SCREW SL 6MM	3.5	6.0	SL
	018-22-C62-33-00	SMART RE Ø3.6X0MM	SMART RE Ø3.6 X 0MM + HEXALOBULAR SCREW RE 0MM	3.6	0.0	RE
	018-22-C62-33-01	SMART RE Ø3.6X1MM	SMART RE Ø3.6 X 1MM + HEXALOBULAR SCREW RE 1MM	3.6	1.0	RE
	018-22-C62-33-02	SMART RE Ø3.6X2MM	SMART RE Ø3.6 X 2MM + HEXALOBULAR SCREW RE 2MM	3.6	2.0	RE
	018-22-C62-33-03	SMART RE Ø3.6X3MM	SMART RE Ø3.6 X 3MM + HEXALOBULAR SCREW RE 3MM	3.6	3.0	RE
	018-22-C62-33-04	SMART RE Ø3.6X4MM	SMART RE Ø3.6 X 4MM + HEXALOBULAR SCREW RE 4MM	3.6	4.0	RE
	018-22-C62-33-05	SMART RE Ø3.6X5MM	SMART RE Ø3.6 X 5MM + HEXALOBULAR SCREW RE 5MM	3.6	5.0	RE
	018-22-C62-33-06	SMART RE Ø3.6X6MM	SMART RE Ø3.6 X 6MM + HEXALOBULAR SCREW RE 6MM	3.6	6.0	RE
	018-22-C62-45-00	SMART RE Ø4.5X0MM	SMART RE Ø4.5 X 0MM + HEXALOBULAR SCREW RE 0MM	4.5	0.0	RE
	018-22-C62-45-01	SMART RE Ø4.5X1MM	SMART RE Ø4.5 X 1MM + HEXALOBULAR SCREW RE 1MM	4.5	1.0	RE
	018-22-C62-45-02	SMART RE Ø4.5X2MM	SMART RE Ø4.5 X 2MM + HEXALOBULAR SCREW RE 2MM	4.5	2.0	RE
	018-22-C62-45-03	SMART RE Ø4.5X3MM	SMART RE Ø4.5 X 3MM + HEXALOBULAR SCREW RE 3MM	4.5	3.0	RE
	018-22-C62-45-04	SMART RE Ø4.5X4MM	SMART RE Ø4.5 X 4MM + HEXALOBULAR SCREW RE 4MM	4.5	4.0	RE
	018-22-C62-45-05	SMART RE Ø4.5X5MM	SMART RE Ø4.5 X 5MM + HEXALOBULAR SCREW RE 5MM	4.5	5.0	RE
	018-22-C62-45-06	SMART RE Ø4.5X6MM	SMART RE Ø4.5 X 6MM + HEXALOBULAR SCREW RE 6MM	4.5	6.0	RE

	018-24-C62-33-01	SMART UF Ø3.6X1MM	SMART UF Ø3.6 X 1MM + HEXALOBULAR SCREW UF 1MM	3.6	1.0	UF
	018-24-C62-33-02	SMART UF Ø3.6X2MM	SMART UF Ø3.6 X 2MM + HEXALOBULAR SCREW UF 2MM	3.6	2.0	UF
	018-24-C62-33-03	SMART UF Ø3.6X3MM	SMART UF Ø3.6 X 3MM + HEXALOBULAR SCREW UF 3MM	3.6	3.0	UF
	018-24-C62-33-04	SMART UF Ø3.6X4MM	SMART UF Ø3.6 X 4MM + HEXALOBULAR SCREW UF 4MM	3.6	4.0	UF
	018-24-C62-33-05	SMART UF Ø3.6X5MM	SMART UF Ø3.6 X 5MM + HEXALOBULAR SCREW UF 5MM	3.6	5.0	UF
	018-24-C62-33-06	SMART UF Ø3.6X6MM	SMART UF Ø3.6 X 6MM + HEXALOBULAR SCREW UF 6MM	3.6	6.0	UF
	018-24-C62-45-01	SMART UF Ø4.5X1MM	SMART UF Ø4.5 X 1MM + HEXALOBULAR SCREW UF 1MM	4.5	1.0	UF
	018-24-C62-45-02	SMART UF Ø4.5X2MM	SMART UF Ø4.5 X 2MM + HEXALOBULAR SCREW UF 2MM	4.5	2.0	UF
	018-24-C62-45-03	SMART UF Ø4.5X3MM	SMART UF Ø4.5 X 3MM + HEXALOBULAR SCREW UF 3MM	4.5	3.0	UF
	018-24-C62-45-04	SMART UF Ø4.5X4MM	SMART UF Ø4.5 X 4MM + HEXALOBULAR SCREW UF 4MM	4.5	4.0	UF
	018-24-C62-45-05	SMART UF Ø4.5X5MM	SMART UF Ø4.5 X 5MM + HEXALOBULAR SCREW UF 5MM	4.5	5.0	UF
	018-24-C62-45-06	SMART UF Ø4.5X6MM	SMART UF Ø4.5 X 6MM + HEXALOBULAR SCREW UF 6MM	4.5	6.0	UF
	018-25-C62-33-01	SMART MF Ø3.6X1MM	SMART MF Ø3.6 X 1MM + HEXALOBULAR SCREW MF 1MM	3.6	1.0	MF
	018-25-C62-33-02	SMART MF Ø3.6X2MM	SMART MF Ø3.6 X 2MM + HEXALOBULAR SCREW MF 2MM	3.6	2.0	MF
	018-25-C62-33-03	SMART MF Ø3.6X3MM	SMART MF Ø3.6 X 3MM + HEXALOBULAR SCREW MF 3MM	3.6	3.0	MF
	018-25-C62-33-04	SMART MF Ø3.6X4MM	SMART MF Ø3.6 X 4MM + HEXALOBULAR SCREW MF 4MM	3.6	4.0	MF
	018-25-C62-33-05	SMART MF Ø3.6X5MM	SMART MF Ø3.6 X 5MM + HEXALOBULAR SCREW MF 5MM	3.6	5.0	MF
	018-25-C62-33-06	SMART MF Ø3.6X6MM	SMART MF Ø3.6 X 6MM + HEXALOBULAR SCREW MF 6MM	3.6	6.0	MF
	018-25-C62-45-01	SMART MF Ø4.5X1MM	SMART MF Ø4.5 X 1MM + HEXALOBULAR SCREW MF 1MM	4.5	1.0	MF
	018-25-C62-45-02	SMART MF Ø4.5X2MM	SMART MF Ø4.5 X 2MM + HEXALOBULAR SCREW MF 2MM	4.5	2.0	MF
	018-25-C62-45-03	SMART MF Ø4.5X3MM	SMART MF Ø4.5 X 3MM + HEXALOBULAR SCREW MF 3MM	4.5	3.0	MF
	018-25-C62-45-04	SMART MF Ø4.5X4MM	SMART MF Ø4.5 X 4MM + HEXALOBULAR SCREW MF 4MM	4.5	4.0	MF
	018-25-C62-45-05	SMART MF Ø4.5X5MM	SMART MF Ø4.5 X 5MM + HEXALOBULAR SCREW MF 5MM	4.5	5.0	MF
	018-25-C62-45-06	SMART MF Ø4.5X6MM	SMART MF Ø4.5 X 6MM + HEXALOBULAR SCREW MF 6MM	4.5	6.0	MF
	018-23-C62-33-00	SMART SH Ø4.5X0MM	SMART SH Ø4.5 X 0MM + HEXALOBULAR SCREW SH 0MM	4.5	0.0	SH
	018-23-C62-33-01	SMART SH Ø4.5X1MM	SMART SH Ø4.5 X 1MM + HEXALOBULAR SCREW SH 1MM	4.5	1.0	SH
	018-23-C62-33-02	SMART SH Ø4.5X2MM	SMART SH Ø4.5 X 2MM + HEXALOBULAR SCREW SH 2MM	4.5	2.0	SH
	018-23-C62-33-03	SMART SH Ø4.5X3MM	SMART SH Ø4.5 X 3MM + HEXALOBULAR SCREW SH 3MM	4.5	3.0	SH
	018-23-C62-33-04	SMART SH Ø4.5X4MM	SMART SH Ø4.5 X 4MM + HEXALOBULAR SCREW SH 4MM	4.5	4.0	SH
	018-23-C62-33-05	SMART SH Ø4.5X5MM	SMART SH Ø4.5 X 5MM + HEXALOBULAR SCREW SH 5MM	4.5	5.0	SH
	018-23-C62-33-06	SMART SH Ø4.5X6MM	SMART SH Ø4.5 X 6MM + HEXALOBULAR SCREW SH 6MM	4.5	6.0	SH

IMPORTANT: The body of Temporary abutment (Titanium Pure Grade 4) and the upper piece based on PEEK are attached and inseparable.

*Each commercial device code consists of the temporary model (titanium body + upper piece in PEEK) and its respective screw (according to the column composition of the commercial model).

The colors of the fastening screw represent the prosthetic interface: PINK interface with SLIM prosthetic components; YELLOW interface with REGULAR prosthetic components; GREEN interface with SHORT prosthetic components; GRAY interface with UNIFIT prosthetic components; BLUE interface with MORSEFIT prosthetic components.

Temporary abutments are compatible with all Plenum Implants, provide the prosthetic interface corresponds to the implant’s internal connection (SL, RE, SH, UF and MF).

Technical Characteristics of Screws: Screws have a thickness of 2 microns of DLC.

HEXAGONAL SCREW: torque is applied to the screw through the hexagonal connection. In this model, the driver transmits torque to the screw with the driver aligned with the screw’s axial axis;

HEXALOBULAR SCREW: torque is applied to the screw through the hexalobular 6-lobe connection. In this model, the driver transmits the torque to the screw with a maximum inclination of 25 degrees between axes. Screw specific to the Smart component model.

2. DEVICE COMPOSITION

Temporary Abutments consist of two (2) parts: 1. a body made from Grade 4 Titanium, according to the specifications provided in ASTM F67 and ISO 5832-2, and PEEK (polyether ether ketone) polymer, according to the specifications contained in ASTM F2026 and ABNT NBR 15723-8; 2. a fastening screw made of titanium alloy (Ti-6Al-4V), according to ASTM F136.

Each commercial code for devices that are part of the Temporary Abutments family consists of the following items: component body, PEEK upper piece, and respective screw, as shown in Table 1 and Figure 1.

FIGURE 1. Schematic model of the Temporary Abutment parts.


It is important to note that the screw is an integral and exclusive part of Temporary Abutment models, and without such an item the device will not fulfill its intended purpose of use.

3. DEVICE PRESENTATION FORM

Temporary abutment [(titanium body + PEEK upper piece) + screw] in the models described in the table above are provided in unique packaging, in sterile conditions, and presented as follow:

- Primary packaging is a polyethylene terephthalate (PET) blister and thermo-sealed with surgical grade paper; containing one (01) unit of the device, sterilized using Gamma Irradiation;
- External packaging is a rigid cardboard box containing the primary packaging and three (03) identification and traceability labels;
- The hospital/health service establishment is responsible for guaranteeing the complete traceability of the device through means of annotation in the patient's medical records, the code and lot of the device used, or identification through one of the three (03) traceability labels that accompany the packaging;
- A printed version of the Instructions for Use (including shipping), can be requested at no additional cost at the telephone number +55 11 3109-9000, or through the electronic address (customerservice@plenum.bio) or through the website (www.plenum.bio/ifu);
- Illustrative images of the primary and external packaging of available commercialized Temporary Abutment models can be found below.

PRIMARY PACKAGING.

EXTERNAL PACKAGING.


ACCESSORIES:

The Temporary Abutment has no accessories.

4. DEVICE RECOMMENDATIONS, PURPOSE, OR INTENDED USE

Temporary Abutments are indicated and used in prosthetic procedures involved in dental implant to prepare and install a temporary dental prosthesis (crown) fastened onto Plenum Dental Implants while awaiting tissue repair after the dental implant's installation or the permanent prosthesis's preparation. These components are primarily intended to act as an intermediate abutment on which a temporary dental prosthesis is fastened, whether in fixed single or multiple restorations involving one or two surgical stages.

Temporary abutments are classified for long-term use and recommended for continuous use lasting more than 30 days. The time for which temporary components are used varies according to the type of rehabilitation in which they will be used. The abutments are conventionally used for more than 30 days, but no established time limit exists for replacing a definitive or temporary component.

NOTE: Temporary abutments must only be used with Plenum Dental Implants (medical devices that are not part of these instructions for use and must be purchased separately, subject to a specific ANVISA Registration). Dental Implants from other brands may result in dimensional incompatibility, improper fitting, micro-movements, and increase the risk of corrosion (electrolytic or galvanic) and fracturing of these devices, compromising the surgical and prosthetic procedure results.

5. PRINCIPLES OF INTENDED DEVICE FUNCTION AND PERFORMANCE

Temporary abutments act as an intermediary element between the dental implant and the prosthesis (crown) to restore the patient's masticatory function and partial aesthetic appearance of the patient's oral cavity. They can be used in immediate or delayed loading procedures (second surgical stage). The prosthetic component is fitted and fastened onto the Plenum Dental Implant by a thread and Morse Taper connection with a hexagonal indexer. The components of this abutment family were designed to attend to different clinical situations, and models are provided in order to allow dentists to select the best option according to: available interocclusal space, the height of existing soft tissue, three-dimensional position of the implant, and type of prosthetic rehabilitation.

6. INSTRUCTIONS FOR USE

- Carefully read the instructions for use, recommendations, and guidelines prior to using Temporary abutments;
- The Temporary Abutments provided in a sterile condition. Verify the device's expiration date before use. Do not use devices if the packaging has been damaged or if the expiration date has expired;
- The techniques for the use and application of the various Temporary abutment models that are part of this family vary according to the preference of the dentist. The dentist is responsible for selecting the surgical protocol, the models and dimensions of the devices, the application technique, and the criteria for monitoring and assessing the surgery results. It is recommended that request diagnostic imaging tests be requested in order to visualize the results obtained anatomically (computed tomography);
- Carefully observe the correct combination and fitting between the Plenum Dental Implant installed in the bone tissue and the respective Temporary abutment platform model selected for use. Both components make use of a Morse Taper connection. Plenum's Dental Implants are not part of these instructions for use, require separate ANVISA registration, and must be purchased separately;
- In order to ensure that Temporary abutments are properly installed and positive results are obtained for surgery, it is essential that the sequence of surgical instruments be used according to surgical protocols recommended by the manufacturer;
- Ensure that the Temporary abutment and surgical instruments are intact and complete prior to starting the surgical procedure;

- Temporary abutments must only be used by properly trained and qualified dental surgeons with proven experience and knowledge of implant dentistry techniques;
- Before the surgical procedure, the dentist must implement rigorous preoperative planning through the use of clinical examinations and images (radiographic and tomographic), defining a surgical plan in advance by selecting the model and size of the device most appropriate to the patient's condition;
- The surgical technique to be used must be indicated by the dentist using their clinical and surgical experience and selecting the most appropriate technique for each patient (individualized procedure). Therefore, a general surgical technique is not recommended for all patients. Surgical protocols must be carried out according to references and the previous experience of the dental surgeon, the Temporary device model that will be used, the technique and installation sequence used, taking conventional and established techniques of implantation and dentistry into consideration;
- Only Temporary abutments and Dental Implants manufactured by Plenum must be used in procedures to guarantee the safety and effectiveness of treatment;
- The device must be handled carefully, particularly during the process of fastening the Temporary abutment onto an osseointegrated implant. The component must be passively seated onto the dental implant platform, thereby avoiding a forced fit, the use of excessive force, or oblique insertion;
- **The maximum recommended torque** for the fastening screw used in this family of components is **20Ncm**.


- **The recommended maximum torque** for the fastening screw used in the temporary prosthesis is **10Ncm**.

- It is highly recommended that imaging tests be used to certify the proper seating of the Temporary abutment onto the dental implant platform;
- Ensure that the Temporary abutment's interocclusal height is appropriate for capturing and preparing the temporary prosthesis. Parts can be customized if necessary;
- Ensure that the parts remaining after a customization are able to withstand the stresses to which they will be exposed;
- Parts must be customized in a clean environment using a high rotation pen, cooling, and irrigation to avoid structural damage to the PEEK components;
- The diameter or thickness of Temporary abutments must not be changed;
- Special care must be taken with the component's hexagonal connection and the driver used to introduce the part, as an unsuitable fit between the two elements may result in damage to the abutment's connection, compromising its integrity and proper fastening to the dental implant;
- Care must also be taken with the surface of the screw and the Temporary Abutment's Morse Taper connection. Both surfaces must be free of particles and solid materials that can prevent perfect fitting and fastening to the implant and temporary dental prosthesis;
- The dentist must carefully assess the repeated introduction and removal of a prosthetic component in the same patient, as this may result in deterioration of the threads, the Morse Taper, and the abutment's connection;
- The dentist must clearly guide the patient, who must understand the importance of follow-up during the postoperative period, in which the stability, integrity, and positioning of the prosthetic component and dental implant used must be verified, as well as the osseointegration and biological fastening of the implant at the appropriate time;
- In order to protect the patient and their surgeon, it is important that the correct and complete traceability of the component used be guaranteed by noting the device code and lot in the patient's records or by attaching to the identification and traceability labels that come with the device to these records.
- Examples of indication: for immediate or delayed loading procedures (second surgical stage) according to the available interocclusal space, height of existing soft tissue, the three-dimensional position of the implant, and type of prosthetic rehabilitation.

IMMEDIATE LOADING

ATTENTION: The figures are merely illustrative. They do not represent the real dimensions and characteristics of the device.

7. PRECAUTIONS

- Dentists must not use the Temporary Abutments until they have carefully read the instructions for use and the information provided on the packaging. All information on the use of instruments and the warnings, precautions, and recommendations included in these instructions for use must be considered. A failure to observe these points may result in complications during procedures;
- Temporary abutments used in SL and RE prosthetic interfaces with reduced diameters (diameters of 3.5 and 3.6 mm) are recommended for use in the rehabilitation of the anterior region (canines and incisors);
- Discard the device if any of the series of packaging that ensures sterility is damaged or opened;
- The instruments' device labels include the date of manufacture, code, device lot number, basic dimensions, ANVISA (Brazilian National Healthy Regulatory Agency) registration number, technician responsible, company information, and the respective symbology that allows the component to be properly identified and traced;
- It is recommended that preoperative examinations be requested according to the knowledge of the patient's condition;
- The dentist must adopt the proper clinical/surgical planning described in the specific implantology literature for immediate and subsequent pre-and postoperative periods;
- Before starting the procedure, the dentist must submit the patient to a rigorous oral antisepsis process to prevent the surgical instrument from coming into contact with contaminating substances;
- Post-operative care after installation of the protective prosthesis and proper oral hygiene is of the utmost importance. The patient must be given instructions regarding the existing risks up until complete osseointegration and biological fastening of the implant and placement of the permanent prosthesis; The patient's ability and willingness to follow instructions from the dentist is one of the most important aspects of the dental implant placement procedure;
- Patient must be warned that a failure to comply with instructions for the post-operative period and period after installation of the temporary dental prosthesis may lead to failure and early loss of the component and implant, requiring additional surgery to remove the dental implant or replace the component;
- If proper osseointegration and biological fastening of the dental implant do not occur at the appropriate time, the dental implant and/or prosthetic component used may fail;
- Temporary abutments must be exclusively handled by properly trained and qualified dentists with proven experience and knowledge of implant dentistry techniques;
- Before using this device, the dentist must understand and plan all stages involved in the surgical procedure and verify whether the model and dimension of the instrument that they intend to use meet requirements for preoperative planning;
- Improper selection of the Temporary abutment that will be used, as well as errors in the indicated use, handling, and installation of the device, may result in damage to the device's structure, which may contribute to failure due to fatigue, fracturing, or loosening of the component;
- If there is an impact, deformation, scratches, cracks, or dents visible in the implant, which may impair the correct and proper functioning of the device, it must be discarded;

- Use of excessive torque during the insertion of a prosthetic component may lead to device failure. Do not exceed 20 Ncm of insertion torque for the fastening screw;
- If a temporary abutment and/or the dental implant is overloaded beyond its functional capacity, excessive bone loss and breakdown of devices and other restoration elements may occur;
- Do not force the introduction of a prosthetic component beyond the depth established by the thread limit, as this may cause deformation in the hexagonal connection of the component and/or the insertion driver and fracture of the device;
- The primary and most likely cause of prosthetic component failures infection: infection, inadequate oral hygiene, mechanical overload, undersized component, initial mechanical instability, use of implant and surgical instrument belonging to other systems, wear in the surgical instrument, lack of previous professional training, contamination of the device, technical error in placement of the component, and a failure to adhere to postoperative instructions on the part of the patient. Dentists must consider the entirety of these aspects during preoperative planning;
- Incorrect handling of instruments with small dimensions inside the patient's mouth may present a risk of aspiration and/or swallowing of these elements;
- The Dentists must rigorously assess the initial mechanical stability of a prosthetic component. If an ideal safe mechanical stability has not been obtained, the surgeon must assess and consider the potential removal and/or replacement of the component;
- Smokers, users of tobacco and drugs, alcoholics, diabetics, or malnourished patients must be warned about the increased incidence of lack of osseointegration, with consequent failure of implants and components;
- Elderly patients, patients with mental health problems, individuals with chemical dependency, or any other difficulty in adhering to postoperative recommendations may represent a higher risk of failure of this device since these patients may ignore instructions for this procedure;
- The success of a procedure for placing a prosthetic component is related to the correct indication, selection, resistance of the dental implant and component, correct surgical technique, as well as the quality of osseointegration, which can be affected by numerous biological, clinical and biomechanical factors, which are external to the quality of the dental implants and components. These factors may limit the service life of devices and consequently result in their failure and early loss.

8. WARNINGS

- Single-use only do not resterilize or reuse.
- Do not use it if the package is opened or damaged or the expiration date has expired.

9. CONTRAINDICATIONS

- Temporary abutments are always used with Plenum Dental Implants (provided separately). Dentists must therefore assess the patient's condition prudently before surgery and use of the components and check for any contraindications reported below. Contraindications may include, but are not limited to:

- A dental implant that is poorly positioned and/or presents insufficient osseointegration to properly support or anchor a dental prosthesis and its components or is subject to strong restrictions on the use of this device;
- Patients with poor bone quality due to osteoporosis or osteopenia, which can cause early failure of the dental implant and component due to lack of osseointegration;
- Patients that present signs of allergy or hypersensitivity to the composition of the raw materials used in the manufacturing of the Components (see "Device Composition");
- Insufficient interocclusal or mesiodistal space for prosthetic rehabilitation;
- Inflammation and/or periodontal infection;
- Recent history of systemic or localized infection;
- Patient with developmental disorders affecting the skeletal system;
- Parafunctional habits;
- Patient with insufficient intermandibular space;
- Occlusion and joint problems that are not treatable;
- Pregnancy and lactation;
- Use in patients with immunological changes and systemic diseases that result in poor healing or a lack of integrity in existing tissue at the implant site;
- Patients that are alcoholics or drug users that cannot correctly follow postoperative instructions;
- Patients with a mental illness or that are unable to follow recommendations and care during the postoperative phase;
- Patients that are unable to maintain proper oral hygiene;
- All contraindications included in known dental and scientific literature for the use of dental implants and their components.

10. MAGNETIC RESONANCE IMAGING

The product contains metal parts that may be affected in an MRI environment. Non-clinical tests and magnetic resonance imaging (MRI) simulations have demonstrated that Plenum products made of metallic material can be used conditionally in MRI (MR Conditional).

A patient with this device can be safely scanned in an MRI system under the following conditions:

Static magnetic field: ≤ 3.0 T

Magnetic field with maximum spatial gradient: 4,000 gauss/cm (40 T/m)

Maximum guideline for MRI systems – Specific Absorption Rate (SAR) for whole-body human exposure: 2.3 W/kg (15 minutes of scanning, normal operating mode);

Under the scanning conditions defined above, the medical device is expected to produce a maximum temperature increase of 0.8 °C after 15 minutes of continuous scanning.

In non-clinical tests, an image artifact originating from the medical device extends approximately 11.3 mm from the device when displayed with an echo gradient pulse sequence in a 3T MRI system.

11. POSSIBLE ADVERSE EFFECTS

Adverse effects are only expected if the device does not follow its instructions for use.

However, all. surgical and prosthetic procedures may not achieve positive results. Every procedure presents risks and the possibility of complications, such as infections, allergic reactions to drugs and anesthetic risks, among others. Possible adverse effects associated with the use of prosthetic components may include, but are not limited to:

- Pain or discomfort due to the presence of the prosthetic component and/or dental implant, or due to the surgical procedure, or subsequent to installation of the dental prosthetic;
- Lack of osseointegration, reabsorption of bone tissue, or even delay in bone healing, which can result in the failure and breakage of the prosthetic component and/or dental implant;
- Loss of primary stability resulting from the bone not healing at the appropriate time, or caused by the improper installation of the device;
- Fracture, displacement, or early or late failure of the component and/or implant due to improper installation technique;
- Failure, displacement, or fracture of the component or dental implant due to excessive loads applied to the devices and/or inadequate bone function;
- Loss of the component or dental implant due to parafunctional habits or poor oral hygiene;
- Allergic reaction to the raw materials used in components;
- Alteration of bone structure due to the imbalance of normal forces/stresses on the bone/implant interface;
- Bone necrosis, insufficient vascularization, bone reabsorption, and insufficient bone formation, which can cause the loss of fastening of the dental implant and consequently the component; these effects are correlated to systemic factors associated with the patient or surgical procedure that is inappropriate, without complying with the instructions for use described for this family of components;
- The patient must be informed and understand that complications are possible in all surgical procedures. The dentist must inform the patient of these complications and the risks of implant surgery.

12. SYMBOLS DESCRIPTION

Table 2 describes the symbols that may be on the device labeling. Please check the physical labeling for the symbols applied to the device. The symbols express information related to the identification and use of the device.

TABLE 5. Symbols used for labeling to express information related to the identification and use of the device.

	Manufacturer
	Country and date of manufacture
	Use-by date
	Batch code
	Catalog number
	Sterilized using Irradiation
	Do not resterilize
	Do not use if package is damaged
	Single sterile barrier system
	Keep away from sunlight
	Keep dry
	Do not reuse
	Consult electronic instructions for use
	Unique device identifier
	Device manufacturing material
	Caution: Federal law restricts this device to sale by or on the order of a physician or dentist. (Notification required by FDA for United States market)
Reprocessing prohibited	Reprocessing prohibited (Brazilian Requirement)
	Restricted Use in Magnetic Resonance Imaging

NOTE: The symbols related to the labeling meet the requirements established in the ISO 15223-1:2021 standard.

13. AUXILIARY INSTRUMENTS

- All surgical instruments developed exclusively for the installation of the Temporary abutments are considered to be auxiliary instruments;
- Plenum instruments are made of S420.10 stainless steel that provides them with high strength and durability, according to requirements specified by NBR 13911; ASTM F899, and ISO 5832-1;
- It is essential that the surgeon exclusively use the introducer drivers manufactured by Plenum to avoid damage to the connection of the components and also to introduction drivers;
- A list of surgical and prosthetic instruments can be found in the

respective printed or electronic communication materials or on Plenum's website (<http://www.plenum.bio>). The use of Plenum instruments must follow the manufacturer's instructions;

- Surgical and Prosthetic instruments are not the main purpose of these instructions for use, that is, they have separate registrations with ANVISA, so they must be purchased separately. Use only instruments manufactured by Plenum, so that there is no incongruence between the fittings and irreversible damage to the components.

14. STERILIZATION

- Temporary abutments are provided in sterile conditions, and sterilized using Gamma Irradiation;
- Before use, check the integrity of the packaging and the expiration date. Do not use devices if the packaging has been damaged with or the expiration date has expired;
- Date of manufacture, expiration date, code, lot number and other information regarding identification and traceability are on the device labeling;
- To maintain sterility, the device must be opened only at the time of use and handled in sterile surgical fields;
- Device packaging must only be opened by qualified personnel trained to perform this procedure.

15. EXPIRATION DATE

- Indicated on device labeling.
- Before use, check the expiration date, do not use the device if it has expired.

16. IDENTIFICATION AND TRACEABILITY

To ensure the traceability of Temporary Abutments and comply with health surveillance requirements, the dentist or someone on their team must keep information about the component used in the patient's medical record.

In cases relating to distributors, the information must be passed on to the manufacturer to complete the traceability cycle. The information necessary for traceability is related to the device used and surgical procedure as follows:

- Name of the dental surgeon;
- Name of hospital/clinic/office;
- Name of manufacturer and distributor
- Date of surgery;
- Code of device used;
- Device lot number;
- Quantity used;
- ANVISA device registration number.

The responsible dentist and their team must use the information in the three

(03) traceability labels supplied with the device, fixing this information in the patient's clinical record to maintain traceability. In addition, this information must be provided to the patient so that they have information about the device used in their procedure.

Labels contain the following information necessary for the traceability of the Temporary Abutment used.

- Manufacturer identification;
- Device name or commercial model;
- Device code;
- Device lot number;
- ANVISA device registration number.

Traceability information is necessary for notification by the health service and/or by the patient to the National Health Surveillance Agency – ANVISA and to the manufacturer, when complaints, technical complaints and/or adverse events occur, for the carrying out of the appropriate investigations.

17. STORAGE, TRANSPORT AND HANDLING CONDITIONS

- The packaging must be intact at the time of receipt, do not use the device if the packaging is damaged or if the expiration date has been exceeded. If the packaging is damaged, the device must be considered non-sterile and discarded;
- The device must be stored and transported in a clean environment, away from sources of heat or humidity, sunlight, and the action of the weather, in order to prevent any damage or changes in the packaging or its physical form that renders its use impossible;
- The device must be stored on shelving or drawers/cabinets to allow it to be easily viewed and handled by the operator. Do not store the device near lamps, so that the information contained in the packaging is not affected;
- The device should only be opened at the time of use in a sterile surgical field;
- The handling of the device should only be done by qualified and trained personnel to perform this procedure.

18. DISPOSAL OF DEVICE

- Discard the device if any of the successive packages that ensure sterility is damaged or opened;
- After using the device or its being removed from the patient, devices that are not subject to analysis or subsequent studies must be disposed of by the health service establishment (hospitals and clinics);
- It is recommended that Temporary Abutments be completely disassembled to ensure it cannot be reused. The health service institution/establishment is solely responsible for disassembling the device to ensure that it is not reused;
- After dismantling, the device must be identified as being unsuitable for use and disposed of in accordance with health service waste management standards (hospitals and clinics) and applicable environmental legislation;
- Temporary Abutments that are not used during clinical procedures must also be discarded following the guidelines above.

19. CUSTOMER SERVICES

- If the device is subject to any adverse event that presents a risk of harm to the patient, generates or has the potential to cause injury or a threat to public health, is out of specifications, or is producing unsatisfactory results, the respective health professional must generate a technical report detailing the specific instance of non-compliance, including claimant information, device code, and lot and patient data, directly notifying Plenum's Customer Service;
- If necessary, send the clean, decontaminated, sterilized, packaged, and correctly identified device, with the description of non-conformities, to the following address: M3 Health Indústria e Comércio de Produtos Médicos, Odontológicos e Correlatos S.A. Avenida Ain Ata, N° 640 - Lot 15 - Quadra B - Jardim Ermida I, Multivias II - Polo Industrial e Logística - Postal Code 13212-213 - Jundiaí - SP - Brazil. Phone: + 55 11 3109-9000; Email: customerservice@plenum.bio.
- ANVISA also makes a channel available to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found on the website of the National Health Regulatory Agency - ANVISA at the electronic address www.anvisa.gov.br, using the NOTIVISA link.
- FDA also provides a channel to report incidents, adverse events, and technical complaints related to the use of devices and services subject to health regulations, which can be found in the FDA Manufacturer and User Facility Device Experience (MAUDE) database.

20. LIMITED WARRANTY

Devices manufactured by Plenum are provided with a factory warranty protecting against manufacturing defects. The occurrence of defects in the device must be reported immediately to the manufacturer, observing the legal deadline and attaching a technical report with the details of the non-compliance with the device. This report must be prepared and signed by a qualified health professional, and relevant ANVISA standards (see the Customer Services section) must be observed in cases involving a risk to the patient's health.

The warranty of devices manufactured by Plenum is conditioned to the strict observation and follow-up of the indications for use, contraindications, precautions, warnings, and other information contained in the instructions for use. The lack of such compliance and follow-up exempts the manufacturer and distributor of the devices from any kind of liability, in particular in the following cases:

- Reprocessing, resterilization, or reuse of the device by the user;
- Improper use or use other than in a surgical setting;
- Device not used according to instructions for use.

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